



Foyomed®



# FOYOMED Medical Instruments

**Foyomed®**

NINGBO FOYOMED MEDICAL INSTRUMENTS CO., LTD.

ADD: 9F, Ningshing Building, No.666 jiangnan Road,  
Hi-Tech Zone, Ningbo, China  
Tel: 86-574-81877253 Fax: 86-574-81877057  
<http://www.foyomed.com> [www.foyocare.com](http://www.foyocare.com)  
E-mail: [sales@foyomed.com](mailto:sales@foyomed.com)

**MFLAB®**

NINGBO MFLAB MEDICAL INSTRUMENTS CO., LTD.

ADD: NO.508, Yindong Road(N), Yinzhou Economic Development Zone,  
315145 Ningbo China  
Tel: 0574-81877852 Fax: 0574-81877057  
<http://www.mflab.com> [www.breathcare.cn](http://www.breathcare.cn)  
E-mail: [products@mflab.com](mailto:products@mflab.com) [breathcare@mflab.com](mailto:breathcare@mflab.com)



FOR YOUR HEALTH  
FOR YOUR SMILE



#### Establishment

FOYOMED is established in the year 2003 as a professional import and export corporation for medical products. It is the headquarters of MFLAB Medical with production area covering 15,500 m<sup>2</sup>. Among it, it is including 100000 grade purification workshops and cleaning room, totally 9,200 m<sup>2</sup>. There are more than 250 workers in the production line.

#### Team Work

FOYOMED has a high efficient and professional team, full of wisdom and passion, with 15 years of working experience. We take the diligent, professional and responsible spirit, committed to provide customers good products and service.

#### Export Reputation

Our FOYOMED products have been exported to over 100 countries worldwide through a network of distributors and subsidiaries, covering North American, South American, European, Asian and African countries.

#### Leading Products

FOYOMED has a wide range of products in business with leading products respiratory products, disposable medical products, medical equipment and laboratory products.

#### Brands

There are five brands for our FOYOMED: Foyomed, Foyocare, FoyoLab, Mflab and Smilecare in distribution.

#### Quality System

FOYOMED is the headquarters of MFLAB Medical which has been passed by GMP, ISO13485, US FDA on-site Audit, Saudi Arabian FDA on-site Audit, Malaysian MOH on-site Audit, TUV CE 0197 certification of more than 40 products, and a number of patents. Our FOYOMED also has been passed by ISO 13485

#### Principle

We FOYOMED adheres to "For Your Smile" principle to strive to achieve excellent quality, responsible price and first-rate service in operation, and cherish every business opportunity.



# Contents

Diagnostics Series	P01~26
Examination Series	P27~31
Weighing Series	P32~35
Rehabilitation Series	P36~44





**FYD1111**

Aneroid Aphygmomanometer  
Euro Type  
Cuff: Nylon D ring 18x14cm  
PVC Bladder/Bulb  
Bag: Vinyl



**FYD 1121**

Arm Type  
Cuff: Nylon D ring 50x14cm  
PVC Bladder/Bulb  
Bag: Nylon



**FYD1119**

Zinc Alloy cover chromeplate  $\Phi$ 68mm  
Rotary handle and gauge  
Fast air release valve  
Fast speed connector  
Cuff: Nylon D ring 48x14cm  
PVC Bladder/Bulb  
Bag: Nylon



**FYD1115**

Palm Type Aneroid Sphygmomanometer  
ABS Plastic cover chromeplate  $\Phi$ 67mm  
Cuff: Nylon D ring 48x14cm  
PVC Bladder/Bulb  
Bag: Nylon



**FYD1116**

ABS Plastic cover chromeplate  $\Phi$ 51 mm  
Cuff: Nylon D ring 48x14cm  
PVC Bladder/Bulb  
Bag: Nylon



**FYD1120**

ABS Plastic case  
Nylon Case  
PVC Color Bulb  
O.D. of Gauge:  $\Phi$ 63mm  
Fluorescent Surface



**FYD1122**

ABS Plastic case  
Nylon Case  
PVC Color Bulb  
P.D. of Gauge:  $\Phi$ 69mm



**FYD1123**

Clock-Type

**Specification:**  
 300mmHg Scale glass tubes  
 Metal case  
 Cotton or nylon cuff  
 With or without D ring is available  
 Artery label  
 Latex bladder with double tubes



**FYD1124**

Desk Type



FYD1131



FYD1132

**FYD1131**

Mercury-Free  
 LCD Display

**FYD1132**

Mercury-Free  
 LED Display



**FYD1134**

Mercury-free  
 LCD Display  
 Available discrepancy: +/- 3mmHg  
 Pressurization: manual by bulb  
 Depressurization:  
 Manual by air release valve  
 Adjustable Height



FYD1135



FYD1136

**FYD1135**

Mercury-Free  
 LCD Display

**FYD1136**

Mercury  
 Sphgmomanpmer

Fully automatic Electronic  
blood pressure Monitor



**FYD1138**

Measure scope: 0 - 300mmHg  
Available with PVC and latex bladder/bulb  
Available with Nylon and Cotton cuff

**FYD1137**

Aneroid Sphygmomanometer  
Clock-type Aneroid Sphygmomanometer with stand.  
Two-piece stand column, easily  
Adjustable from 80 to 125cm, 5pcs castors with brakes,  
individual packed in foam box.



**Wrist Type**



**FYD1151**

LCD Display 35x48mm  
7 sets including blood pressure,  
pulse rate, date, time  
Measurement Range: 0 - 300mmHg  
40 - 200Pulse/min(pulse)  
Battery: AAA\*2



**FYD1154**

Large LCD display with bold type  
90 sets of memories



**FYD1155**

Large LCD Display  
Five language voice adjustable(E. F. R.G.S)  
Voice volume adjustable  
Two users with 90 sets of memory each  
Automatic power off after 3 minutes  
of no operation



## Arm Type

### FYD1150

LCD Display 80x23mm  
 21 sets including blood pressure,  
 pulse rate, date, time  
 Measurement Range: 0 - 300mmHg  
 40 - 200 Pulse/min(pulse)  
 Battery: AAA\*4



### FYD1153

Extra large LDC display and back light  
 Five language voice adjustable(E. F. R.G.S)  
 Voice volume adjustable  
 Two users with 90 sets of memory each  
 Average for the latest 3 reading



### FYD1152

Backlighting in 3 different colors  
 Intelligent inflation  
 Universal Cuff 22 -42cm  
 LCD Backlight display  
 90 sets of memory  
 Battery: AAA\*4



## Single Head Stethoscope



Oxidized or Anodized Aluminum chest  
 Brass Chrome Plated Binaural  
 Plastic eartips  
 All kinds of color available for the tube such as:  
 Black, Red, Green, Yellow, Blue, Purple, Pink,  
 Grey, White, etc.



### FYD1210A

Single head stethoscope  
 Colored chest piece



### FYD1211

Single head stethoscope  
 Silver chest piece



### FYD1212

Single head stethoscope  
 Colored chest piece & Binural



### FYD1213

Single head stethoscope  
 with plastic ring

**Stethoscope**

1. Material of chest piece: Zinc Alloy, aluminum Alloy, Stainless steel
2. Material of earloop: Copper, Aluminum, Stainless steel.
3. Material of earplug: PVC, ABS
4. Material of tube: PVC.
5. Any color of tube is available.



**FYD1215**

single head stethoscope  
Zinc alloy chest piece  
with chrome plated silver finish  
Deluxe Special



**FYD1214**

Pediatric  
single head stethoscope



**FYD1216**

Bowles stethoscope  
Chrome plated chest piece



**FYD1217**

Disposable Single head stethoscope  
Plastic



**Dual Head Stethoscope**  
Oxided or Anodized aluminum chest piece  
Brass chrome plated binaural  
Plastic or soft eartips  
All kinds of color available for the such as:  
Black, Red, Green, Yellow,  
Blue, Purple, Pink,  
Grey, White,  
etc.



**FYD1210B**

Dual head stethoscope  
Colored chest piece



**FYD1221**

Dual head stethoscope  
Colored chest piece & Binaural

**FYD1220**

Sliver chest piece



**FYD1222**

Pediatric  
Dual Head Stethoscope



**FYD1223**

Teaching use  
dual head stethoscope



**FYD1224**

Neonatal  
Dual Head Stethoscope



**FYD1225**

Dual head stethoscope  
with non-chill ring with soft eartips



**FYD1226**

Dual head stethoscope with non-chill ring  
with soft eartips



**Stainless steel stethoscope**

Stainless steel chest piece  
Stainless steel Binaural  
Plastic & Soft eartips  
With non-chill rings  
Spare diaphragm included  
All kinds of color available for the such as:  
Black, Red, Green, Yellow, Blue, Purple,  
Pink, Grey, White, etc.



**FYD1232**

Stainless steel stethoscope

**FYD1233**

Stainless steel stethoscope  
cardiology type



**FYD1230**

Stainless steel stethoscope  
Adult type ingle-lumen tubing



**FYD1231**

Stainless steel stethoscope  
Child type



**FYD1242**

Master Cardiology Stethoscope

Specification:

Application: Cardiology/excellent performance.

Product weight/length: 215g/71cm

Chest piece material:stainless steel,

surface treatment: manual polish

Diaphragm:Φ44mm, Type: floating diaphragm

Binaural material: stainless steel, Φ6mm

Eartip material: silicone, structure: Soft sealing structure

Accessory box: 2 sets silicone ear tips

Package Content: 1pc/box, 12pcs/ctn, 580x400x245mm



Sprague Rappaport Stethoscope



**FYD1241**

Sprague Rappaport type chest piece  
Brass chrome plated binaural  
Plastic or soft eartips  
All kinds of color available for the such as:  
Black, Red, Green, Yellow, Blue, Purple,  
Pink, Grey, White, etc.  
2sets of earplug, 1set of mebrane  
1pc of large chest piece  
1pc of middle chest piece  
1pc of small chest piece



**FYD1243**

Fetal Stethoscope  
Metal Type



----- Accessories -----

**FYD1250**



**FYD1253**

Earplug



**FYD1101**

**Buck Neurological Hammer**  
Hammerhead: zinc alloy with rubber plug  
Handle: Nylon  
With needle: Zinc alloy, with brush  
Length of hammer: 22.5cm  
Height of hammer head: 8.5cm



**FYD1102**

**Gear Wheel Hammer**  
Hammerhead: Zinc Alloy  
Handle: Zinc Alloy



**FYD1103**

**Taylor Hammer**  
Length of Handle: 18cm  
Size of Runner Plug: 5 x 4.3 x 1 cm  
Thickness of Rubber Plug: 1.5cm

**FYD1105**

**Gristletly Hammer**  
Rubber Plug  
Hammerhead: copper  
Handle: steel  
Length of hammer: 24cm  
Diameter of hammerhead: 4.5cm



**FYD1104**

**Babinski(Adult) with Telescoping Handle**  
Hammerhead: Zinc alloy with rubber plug  
Handle: Zinc Alloy  
Length of handle: 21cm  
Diameter of hammerhead: 4.5cm  
Thickness: 1cm



**FYD1106**

**Colorflex reflex hammer**  
Hammerhead: rubber with copper ring  
Handle: Nylon  
Length of hammer: 22.5cm  
Height of hammerhead: 7.5cm





**FYD1201**

Diagnostic Pen Light

**FYD1202**

Diagnostic Pen Light



**FYC1509**

Specification:

- 1. One chargeable Otoscope & Ophthalmoscope
- 2. Single Charger
- 3. Gift box

HS-OTP10C-2

Specification:

- 1. One chargeable Otoscope & Ophthalmoscope
- 2. Double Charger
- 3. Gift box

**FYC1510**

Integrated Diagnostic System

- Fiber Optic Otoscope
- Direct Ophthalmoscope
- Tip Dispenser(with disposable tips)
- Infrared thermometer
- Insufflation Bulb
- Sphygmomanometer
- Customize your wall set to match your individual needs



**FYC1501**

ENT Kit

- Ophthalmoscope 1
- Otoscope carrier 1
- Ear specula(3,4,5,6mm diam) 4
- Nasal apeculum(8mm diam) 1
- Dilating Nasal speculum 1
- Laryngeal mirror 3
- Tongue depressor plastic 1
- Battery handle 1
- Spare Bulbs(2large, 4small) 6



**FYC1502**

Professional Otoscope

- Features:
- Fiberoptics for cool, obstruction free light transmission
- Otoscope head made from polyimid - an impact resistant space age thermoplastic Swivelling viewing window with 2.5x magnification
- "AA" polycarbonate battery handle interchanges pocket clip with intergral on/off slide switch for ambidextrous use
- Technical Data
- Magnification: 2.5x
- Distal dia.: 2.5mm, 4mm
- Bulb: 2.5V 1.6W halogen lamp
- Batteries: two "AA" Batteries or two "AA" NiMH Rechargeable Battery



**FYC1503**

Otoscope

Deluxe packing:  
Penlight, Part of otoscope, batteries



**FYC1504**

Reflex Hammer Gifts Set

Buck Neurological Hammer  
Babinski with telescoping handle  
Taylor Hammer



**FYC1506**

Buck neurological hammer

Gear wheel hammer  
Taylor Hammer  
Gristletly hammer



**FYC1508**

Fiber laryngeal mitrrior

Power supply: 1.5V X 2 "C" Type  
Bulb: 2.5V/0.35A  
Blades: 92, 114, 134mm



**FYD1465**

Non-contact Infrared Forehead thermometer  
 Multiple USE ( 3 in 1 ) :  
 a. forehead temperature  
 b. surface of object temperature  
 c. room temperature  
 Display Range & Accuracy:  
 a. forehead temperature: 32°C-42.9°C , ±0.3°C;  
 b. surface of object temperature: 0°C-60°C , ±1°C;  
 c. room temperature: 0°C-40°C , ±0.5°C  
 Measuring time: 1 second  
 Reading store: 32 data memory  
 Battery: 3.0V (AAA\*2)  
 Automatic switch-of  
 Tricolor illuminated display (fever alarm)  
 Voice broadcast (speaking)



**FYD1467**

Specification:  
 waterproof ( Test Method: located at least 15 cm below surface of water for 30 min ).  
 Flexible tip.  
 10 seconds.  
 Beeper & fever alarm  
 Last readout recall  
 automatic switch-off  
 low battery indication  
 accuracy ±0.1°C(35.5°C-42.0°C).  
 Battery:1.5V(LR/SR41).

**FYD1464**

Infrared Thermometer  
 Range of displayed temperature: 32.0°C - 42.9°C(89.6°F - 109.2°F)  
 Ambient temperature operating range: 10°C-40°C(50°F - 104°F)  
 Storage temperature range: -25°C -55°C(-13°F -131°F)  
 Digital increment: 0.1°C/0.1°F  
 Maximum permissible error:  
 ±0.2°C/0.4°F 35.5 - 42.0°C(95.9°F -107.8°F)  
 ±0.3°C/0.5°F 32.0 - 35.4°C(89.6°F -95.7°F)  
 ±0.3°C/0.5°F 42.1 - 42.9°C(107.8°F -109.8°F)



**FYD1425**

Specification:  
 Transparent.  
 60 seconds.  
 Beeper & fever alarm  
 Last readout recall  
 automatic switch-off  
 low battery indication  
 accuracy ±0.1°C (35.5°C-42.0°C) .  
 Battery:1.5V(LR/SR41).



**FYD1460**

Infrared Ear Thermometer  
 3-Digital  
 Scale: °C or °F  
 Accuracy: ±0.2°C(±0.4 °F)  
 Min Scale: 0.1 °C  
 Memory function: 10 group  
 Measuring times: 1 second  
 Beeper function  
 Auto shut-off



**FYD1466**

waterproof ( Test Method: located at least 15 cm below surface of water for 30 min )  
 Backlight  
 Jumbo LCD: 33.50X16.0 mm  
 Beeper & fever alarm  
 Last readout recall  
 automatic switch-of  
 low battery indication  
 accuracy ±0.1°C(35.5°C-42.0°C)  
 Battery:3.0V(CR2032)



**FYD1437**

Specification:  
10 seconds  
Beeper & fever alarm  
Last readout recall  
automatic switch-of  
low battery indication  
accuracy  $\pm 0.1^{\circ}\text{C}$  ( $35.5^{\circ}\text{C}$ - $42.0^{\circ}\text{C}$ ).  
Battery: 1.5V(LR/SR41).



**FYD1441**

waterproof  
flexible tip.  
10 seconds.  
Beeper & fever alarm  
Last readout recall  
automatic switch-of  
low battery indication  
accuracy  $\pm 0.1^{\circ}\text{C}$  ( $35.5^{\circ}\text{C}$ - $42.0^{\circ}\text{C}$ ).  
Battery: 1.5V(LR/SR41).



**FYD1422**

Large LCD display  
Display range:  $32.0^{\circ}\text{C} \sim 42.0^{\circ}\text{C}$   
 $89.6^{\circ}\text{F} \sim 107.6^{\circ}\text{F}$   
Accuracy:  $\pm 0.1^{\circ}\text{C}$ ,  $\pm 0.2^{\circ}\text{C}$   
Min scale:  $0.1^{\circ}\text{C}/0.1^{\circ}\text{F}$   
Response: 60seconds  
Size:  $121 \times 17 \times 10\text{mm}$   
LCD:  $20 \times 7.5\text{mm}$   
N.W.: 9.5g



**FYD1434**

Flexible Probe  
Water-proof  
Large LCD Display  
Display Range:  $32^{\circ}\text{C} - 42^{\circ}\text{C}$   
 $89.6^{\circ}\text{F} - 107.6^{\circ}\text{F}$   
Accuracy:  $\pm 0.1^{\circ}\text{C}$ ,  $\pm 0.2^{\circ}\text{F}$   
Min Scale:  $0.1^{\circ}\text{C}/0.1^{\circ}\text{F}$   
Response: 60seconds  
Size:  $126 \times 18 \times 10\text{mm}$   
LCD:  $20 \times 7\text{mm}$   
N.W.: 11.5g

Thermometer

Clinical Thermometer

Glass  
Scale:  $^{\circ}\text{C}$ ,  $^{\circ}\text{C}$  &  $^{\circ}\text{F}$   
Accuracy:  $\pm 0.1^{\circ}\text{C}$  ( $\pm 0.2^{\circ}\text{F}$ )  
White back, Yellow back or Blue Back



**FYD1311**

Oral Use



**FYD1312**

Armpit Use, Large Size



**FYD1313**

Armpit use, Middle Size



**FYD1320** Air Thermometer

**FYD1330** Bath Thermometer



**FYD1451**

Bay Pacifier  
Digital Thermometer(Waterproof)  
Display range:  $32.0 - 42.0^{\circ}\text{C}$  /  $90.0 - 107.6^{\circ}\text{F}$   
Accuracy:  $\pm 0.1^{\circ}\text{C}/\pm 0.2^{\circ}\text{F}$   
Min. Scale: 0.1  
Memory: Last measurement reading  
Measurement time(Reference only,  
it differs from people(to people): 2min  
Beeper function  
Auto shut-off  
Battery: 1.5V button battery(LR/SR-41)  
Size:  $62 \times 46 \times 50\text{mm}$   
LCD:  $15 \times 6\text{mm}$   
N.W.: 13G  
Memory: last measuring reading

**FYD1501**

Alcohol Tester  
Advanced NM Hot-wire alcohol sensor  
Quick Response  
SMD Assembling, Stable performance  
Smart MCU control  
Direct testing process LCD Indication  
Portable design  
No touch health Design  
Audio and visible warning beyond preset level  
Battery saving design, low voltage indication



**FYD1401**

Height measuring Rod

Height range to be measured: 0-2000mm  
 Min value of height per division: 1mm  
 Measurement: 35.5×19×19cm /20pc  
 G.W.: 4kg N.W.: 3.6kg



**FYC1701**

Luxurious Eyesight Lamp box

This is served as an apparatus for enghtsight examination in hospitals In practical in hospitals Soft and comfortable lighting levels and being made of aluminum materials, it has a very handsome appearance There are three styles available-desktop style, wall style and hanging style  
 Main Technical data:  
 Power Voltage: AC 220V±10%  
 Frequency: 50Hz  
 Power: 60W



**FYC1801**

Operation Reflector Lamp

Stainless haulm,  
 Hose flexible lift type, unbending rise or drop  
 Voltage: AC 220V  
 Frequency: 50Hz  
 Power: 60W



**FYW1101**

Electronic Infant Scale  
 Max. Weight: 20kg



**FYW1102**

Baby Scale  
 Max. Weighing: 20kg



**FYW1103**

Electronic Infant Scale  
 Max. Weight: 20kg  
 Min. Value per division: 50g



**FYW1104**

Max. Weighing: 25kg



**FYW1105**

Baby Scale  
 Tray Size: 56x28x9.3cm  
 Body Size: 36x29.5x7.5cm  
 Net Size: 56x31.5x13.5cm  
 Max. Weight: 30kg  
 Graduation: 5g  
 LCD Display, Auto Power off,  
 Low battery indicator, 1 loadcell  
 Precision: ±1%  
 Displayed In two units: KG/LB/OZ  
 Batteries: 4X AA(not included)  
 LCD Size: 7 x 2.8cm



**FYW1301**

Bathroom scale



**FYW1302**

Bathroom scale



**FYW1303**

Bathroom scale



**FYW1304**

Bathroom scale



**FYW1305**

Bathroom scale



**FYW1306**

Bathroom scale



**FYW1201**

Health Scale  
Max. Weighing: 120kg



**FYW1203**

Health Scale  
Max. Weighing: 200kg

**FYW1204**

Electronic Body Scale  
Max. Weight: 200kg  
Division: 100g  
Display: LED Power: AC/DC  
Height range to be measured: 800 - 2100mm  
Min value of height per division: 5mm  
G.W.: 15KG N.W.: 13kg  
Measurement: 100x32x26cm





FYR1102












FYR1103



FYR1104

**MAIN PARAMETERS**

MEASUREMENTS(cm) MODEL NO.							 (Kg)	 DIM	 NW/GW
FYR1102	64	46	58	20	49	87	100	94X23X89	16.7 19
FYR1103	64	46	58	17	48	88	100	92X28X89	15.3 16.5
FYR1104	65	46	58	20	51	119	100	98X35X89	27.7 29.7



FYR1106



FYR1110



FYR1114



FY1111



FYR1107



FYR1139

MAIN PARAMETERS

MEASUREMENTS(cm) MODEL NO.							 (Kg)	 DIM	 NW/GW	
FYR1106	65	46	62	20	51	88	100	80x28x89	19.3	22
FYR1112	52	35	58	15	48	90	75	85x33x73	13	15.5
FYR1110	64	44	58	20	55	129	100	107x33x78	25.6	29.3

MAIN PARAMETERS

MEASUREMENTS(cm) MODEL NO.							 (Kg)	 DIM	 NW/GW	
FYR1139	38	32	15	10	45	84	100	39x29x70	6	7.5
FYR1114	61	46	51	15	49	92	100	82x28x76	11.7	13.7
FYR1111	61	46	56	17.5	49	58	100	90x27x91	11.8	14.1





FYR1113



FYR1115

MAIN PARAMETERS

MEASUREMENTS(cm) MODEL NO.									
FYR1113	52	38-44	38	15	49	96-114	75	77x40x95	21.4 24.3
FYR1115	52	38-44	40	15	48	97-116	72	84x40x95	18.3 21



FYR1116



FYR1117

MAIN PARAMETERS

MEASUREMENTS(cm) MODEL NO.																		
FYR1116	63	46	32	20	50	91	20	1-7	28AH	200Wx2	4°	12°	100	80x40x76	59.7	66.2		
FYR1117	64	46	32	20	55	97	20	1-7	28AH	200Wx2	4°	12°	100	78x68x69	62.8	70.8		



FYR1201



FYR1202



FYR1203



FYR1204



FYR1205

MAIN PARAMETERS

MODEL NO.								
FYR1201	47	46	78-96	53X17X81	4.8	2	5.7	
FYR1202(S)	43	38	65-75	76x68x35	14.9	10	16.9	
FYR1202(M)	44	38	75-85	86x68x35	15.2	10	17.5	
FYR1202(L)	45	38	85-95	96x68x35	16	10	18.7	
FYR1203	47	47	78-96	53X17X81	5.1	2	6.2	
FYR1204	50	44	78-96	58x18x82	5.5	2	6.7	
FYR1205	39	62	51-61	52.5x22x63	2.9	2	3.9	

WALKER

Commode Chair-without wheel



FYR1301



FYR1302

MAIN PARAMETERS

Item Code									PCS CN	G.W (kg)
FYR1301	56	46	41-51	74-84	40	33	6	60X32X48	2	13.4
FYR1302	41	38	40		40		1.6	75X30X40	6	10.7

STICK



FYR1401



FYR1402



FYR1403



FYR1404



FYR1405

MAIN PARAMETERS

MODEL Number	HEIGHT (cm)	N.W. (Kg)		PCS CN	G.W. (Kg)
FYR1401(L)	128-148	0.9	132x33x29	20	19
FYR1401(M)	116-136	0.85	119x33x29	20	18
FYR1401(S)	103-123	0.8	103x33x29	20	17
FYR1405	82-92	0.3	64x29x16	20	7
FYR1402	77.5-100	0.3	79x29x18	20	7
FYR1404	120	0.3	35x20x26	40	13
FYR1403	77.5-100	0.35	79x34x19	20	8

FYR1408



FYR1406



FYR1407



FYR1409



FYR1410



FYR1411



FYR1412

MAIN PARAMETERS

MODEL Number	HEIGHT (cm)	N.W. (Kg)		PCS CN	G.W. (Kg)
FYR1406	89	0.3	94x28x15	20	7
FYR1410	72.5-95	0.9	79x34x43	10	10
FYR1411	72.5-95	0.9	60x28x47	10	10
FYR1409	65-87.5	0.8	79x34x43	10	9
FYR1408	68.5-91	0.3	65x28x15	20	7
FYR1407	68.5-91	0.3	65x28x15	20	7
FYR1412	80-102.5	0.9	60x28x47	10	10



**FYR1601**

Cervical Collar  
Macromolecule Foam Materials, light



**FYR1502**

Tractor For Neck Ache  
Designer for treating and preventing neck fractured is location of vertebra. Protrusion of the intervertebral disc and excessive stretch and draw. It can support the patient's head efficiently depending on appropriate air pressure, so that it can release the neck and shoulder burden. Also it has the magetic infrared treatment function to dispel the pain.



**FYE1302 L**

Head Immobilize  
Unfolded Size(LxWxH)  
40x26x16cm  
Self Weight: 1Kg  
Packing: 8pcs/54x45x63cm  
G.W.: 11Kgs



**FYE1302 C**

Head Immobilize  
Unfolded Size(LxWxH)  
26x34x14cm  
Self Weight: 0.5Kg  
Packing: 10pcs/58x38x56cm  
G.W.: 8Kgs



**FYE1303 A**

Cervical Collar  
Unfolded Size(LxWxH): 59x19x0.5cm  
Self Weight: 0.2kg  
Packing: 40pcs /59x43x 27cm  
G.W.: 7kgs



**FYE1303 B**

Cervical Collar  
Unfolded Size(LxWxH): 54x13x0.5cm  
Self Weight: 0.1kg  
Packing: 40pcs /61x45x 39cm  
G.W.: 6kgs



**FYR1603**

Cervical Collar  
Size: S.M.L





**FYU1501**

**Mattress Characteristics:**

Mattress with adjustable pump  
 Max. User Weight : 160 Kgs.  
 Mattress Dimensions : 185 x 80x9.5 cm.  
 Materials :Nylon PVC cell  
 Cycle Time : 6 Mins.  
 Air Output : 6-7 L/ Min.  
 Pressure Range : 70~ 130 mmHg  
 Electrical Power Supply : 220-240 Volts.  
 Protection Against Class II Electric Shocks  
 Warranty : 12 months



**FYU1502**

**Mattress Characteristics:**

Mattress with adjustable pump  
 Max. User Weight : 130 Kgs.  
 Mattress Dimensions : 200 x 90 cm.  
 Materials : PVC 0.3mm.  
 Cycle Time : 12 Mins.  
 Air Output : 5-7 L/ Min.  
 Pressure Range : 50~ 110 mmHg  
 Protection Against Class II Electric Shocks  
 Warranty : 12 months



**FYU1503**

Work Mode	Continuous Wave
Voltage	AC220V/110V 50/60HZ
Power	≤10W
Pressure Range	≥16kPA
Air output	≥4.5L/min
Cycle Time	5 - 6minutes
Dimension	1950x900mm
Height after charge	7cm



**FYU1504**

Work Mode	Continuous Wave
Voltage	AC220V/110V 50/60HZ
Power	≤10W
Pressure Range	10 - 18Kpa
Air output	≥4L/min
Cycle Time	6 - 8minutes
Dimension	1950x850mm
Height after charge	9cm



**FYR1501**

Professional Dental Oral Irrigator

Specification:  
 Rated Voltage: 230V AC~50Hz  
 Power consumption: 40W  
 Water Pressure: 90~200Psi  
 Pump frequency: 1200imp/min  
 Max operating time: 10 minutes  
 Water tank: 500ml  
 Dimension: 180mmx165mmx140mm  
 Recommended by dental professional for adults and children  
 For prevention of plaque gingivitis and many dental diseases.  
 Package content: 1pc/box, 8pcs/carton, 65x39x23.5cm.  
 G.W./N.W.: 9.1/8.3kgs.

**FYR1713**

Professional Dental Oral Irrigator

Specification:  
 Rated Voltage: 230V AC~50Hz  
 Power consumption: 40W  
 Water Pressure: 90~200Psi  
 Pump frequency: 1200imp/min  
 Max operating time: 10 minutes  
 Water tank: 500ml  
 Dimension: 180mmx165mmx140mm  
 Recommended by dental professional for adults and children  
 For prevention of plaque gingivitis and many dental diseases.  
 Package content: 1pc/box, 8pcs/carton, 65x39x23.5cm.  
 G.W./N.W.: 9.1/8.3kgs.



**FYR2101**

Oxygen Absorber

Service Medium: O<sub>2</sub>  
 Metering Range(L/min): 1-10  
 Accuracy: grade 4  
 Entrance Pressure(MPa): 12Mpa  
 Exit Pressure(MPa): 0.2-0.3MPa  
 Entrance Connection Screw thread: G5/8-14-RH  
 Output Connection Screw thread: 8mm





**FYR1705**

Ultrasonic Nebulizer

Specification:  
 Power: 220 - 240VAC 50/60Hz  
 Power consumption: Approx. 40VA  
 Max. Nebulizing rate:  $\geq 2.0$ ml/min  
 Particle size:  $\approx 5\mu$ m  
 Medication cup capacity: 12ml  
 Noise Level:  $\leq 55$ dB  
 Carton Size: 520x500x225mm  
 N.W.: 10.6kgs G.W.: 11.6kgs

**FYR1701**

Compressor Nebulizer

Specification:  
 Rated Voltage: AC 230V, 50Hz  
 Voltage Adjustable: 230V/120V(Optional)  
 Extreme Pressure:  $> 250$ Kpa  
 Nebulizing Pressure: 100 - 190Kpa  
 Noise Level:  $\leq 60$ dB  
 Max Capacity of nebulizer kit: 10ml  
 Nebulizing Rate:  $\geq 0.2$ ml/min  
 MMAD: 4.3 $\mu$ m  
 FPD:  $\geq 70\%$   
 Continuous working time:  $\geq 4$  hours  
 Carton size: 390x320x440mm  
 N.W.: 9.8kgs G.W.: 10.5kgs



**FYR1710**

Compressor Nebulizer

Specification:  
 Rated Voltage: AC 230V, 50Hz  
 Voltage Adjustable: 230V/120V(Optional)  
 Extreme Pressure:  $> 205$ Kpa  
 Nebulizer Pressure: 100 - 190Kpa  
 Noise Level:  $\leq 60$ dB  
 Max Capacity of nebulizer kit: 10ml  
 Nebulizing Rate:  $\geq 0.2$ ml/min  
 MMAD: 4.5 $\mu$ m  
 FPD:  $\geq 70\%$   
 Continuous working time:  $\geq 4$  hours  
 Carton size: 550x410x245mm  
 N.W.: 8.3kgs G.W.: 9.1kgs



**FYR1711**

Compressor Nebulizer

Specification:  
 Rated Voltage: AC 230V, 50Hz  
 Voltage Adjustable: 230V/120V(Optional)  
 Extreme Pressure:  $> 205$ Kpa  
 Nebulizer Pressure: 100 - 190Kpa  
 Noise Level:  $\leq 55$ dB  
 Max Capacity of nebulizer kit: 10ml  
 Nebulizing Rate:  $\geq 0.2$ ml/min  
 MMAD: 4.2 $\mu$ m  
 FPD:  $\geq 75\%$   
 Continuous working time:  $\geq 4$  hours  
 Carton size: 305x230x275mm  
 N.W.: 6.0kgs G.W.: 7.0kgs



**FYR1704**

Compressor Nebulizer

Specification:  
 INPUT: 100 - 240V, 50 - 60Hz  
 Rated Voltage: DC 12V 1A  
 Extreme Pressure:  $> 125$ Kpa  
 Noise Level:  $\leq 55$ dB  
 Max Capacity of nebulizer kit: 10ml  
 Nebulizing Rate:  $\geq 0.2$ ml/min  
 MMAD: 3.5 $\mu$ m  
 FPD:  $\geq 85\%$   
 Continuous working time: 20 minutes on, 40 minutes off  
 Carton size: 470x390x270mm  
 N.G.: 9.5kgs G.W.: 10.2kgs

**FYR1712**

Compressor Nebulizer

Specification:  
 Rated Voltage: AC 230V, 50Hz  
 Voltage Adjustable: 230V/120V(Optional)  
 Extreme Pressure:  $> 205$ Kpa  
 Nebulizer Pressure: 100 - 190Kpa  
 Noise Level:  $\leq 60$ dB  
 Max Capacity of nebulizer kit: 10ml  
 Nebulizing Rate:  $\geq 0.2$ ml/min  
 MMAD: 4.2 $\mu$ m  
 FPD:  $\geq 75\%$   
 Continuous working time: 20 minutes on, 40 minutes off  
 Carton size: 400x290x290mm  
 N.W.: 7.0kgs G.W.: 7.9kgs





**FYD151007**

PE Material, 1000cc  
Size: 32 x 21 x 10cm  
60g



**FYD151008**

PP Material  
Size: 36.5 x 28 x 10.6cm  
180g



**FYD151009**

PE Material  
Size: 34 x 26.5 x 9cm  
300g



**FYD152001**



**FYD152002**



**FYD152003**



**FYD152004**



**FYD151101**

News Papper  
Size: 25.5 x 11.5 x 11cm



**FYD151102**

News Papper  
Size: 25 x 11 x 10.8cm



**FYD151103**

PP Material  
Size: 25 x 10 x 4.9cm  
25g



**FYD151011**



**FYD151012**



**FYD151013**



**FYD151014**



**FYD151004**



**FYD151005**



**FYD151006**



**FYD151001**



**FYD151002**



**FYD151003**