

Foyomed[®]



For your smile

Disposable Medical Instruments

Ningbo Foyomed Medical Instruments



Foyomed[®]

NINGBO FOYOMED MEDICAL INSTRUMENTS CO., LTD.

ADD: 9F, Ningshing Building, No.666 jiangnan Road,
Hi-Tech Zone, Ningbo, China

Tel: 86-574-81877253 Fax: 86-574-81877057

<http://www.foyomed.com> www.foyocare.com

E-mail: sales@foyomed.com

MFLAB[®]

NINGBO MFLAB MEDICAL INSTRUMENTS CO., LTD.

ADD: NO.508, Yindong Road(N), Yinzhou Economic Development Zone,
315145 Ningbo China

Tel: 0574-81877852 Fax: 0574-81877057

<http://www.mflab.com> www.breathcare.cn

E-mail: products@mflab.com breathcare@mflab.com

ABOUT US

OUR COMPANY



Establishment

FOYOMED is established in the year 2003 as a professional import and export corporation for medical products. It is the headquarters of MFLAB Medical with production area covering 15,500 m². Among it, it is including 100,000 grade purification workshops and cleaning room, totally 9,200 m². There are more than 350 workers in the production line.

Team Work

FOYOMED has a high efficient and professional team, full of wisdom and passion, with more than 20 years of working experience. We take the diligent, professional and responsible spirit, committed to provide customers good products

Export Reputation

Our FOYOMED products have been exported to over 100 countries worldwide through a network of distributors and subsidiaries, covering North American, South American, European, Asian, Oceanian and African countries.

Leading Products

FOYOMED has a wide range of products in business with leading products respiratory products, disposable medical products, medical equipment and laboratory products.

Brands

There are five brands for our FOYOMED: Foyomed, Fyocare, FoyoLab, MFlab and Smilecare in distribution.

Quality System

FOYOMED is the headquarter of MFLAB Medical which audited by TUV Rheinland every year, more than 30 products has got CE certificates issued by TUV Rheinland. MFLAB also passed by GMP, ISO13485, US FDA on-site Audit, Saudi Arabian FDA on-site Audit, Malaysian MOH on-site Audit, Korea MOH on-site Audit.

Principle

We FOYOMED adheres to " For Your Smile" principle to strive to achieve excellent quality, responsible price and first-rate service in operation, and cherish every business opportunity.

Contents

Breathcare Anesthesia	P01~23
Urology	P24~35
Hypodermic	P36~52
Examination	P53~58
Wound Dressing	P59~73
Medical Non-woven & Accessories	P74~83
Others	P84~92



BREATHCARE ANESTHESIA

Spacer for Aerosol

The spacer is mainly used for children and patients with poor coordination of inhalation and hand movements, or patients with obvious local side effects in the throat after using MDI. Especially for inhaled corticosteroids, the use of a spacer significantly reduces the local and potentially systemic side effects of inhaled corticosteroids.



LB2122

Chamber Volume: 170ml
Square or Round Mouthpiece is available

Cat. NO.	Size	Applicable People	Color
LB212201	S	<18 Months	Orange, Yellow, Dark blue, Green, Whithe
LB212202	M	1-5 Years	
LB212203	L	5+ Years	
LB212204	Mouthpiece	5+ Years	

Patent Number : CN 219049891U



LB2123

Chamber Volume: 200ml, 250ml, 350ml, 500ml
Square or Round Mouthpiece is available

Cat. NO.	Size	Applicable People	Color
LB212301	S	<18 Months	Orange, Yellow, Dark blue, Green, Whithe
LB212302	M	1-5 Years	
LB212303	L	5+ Years	
LB212304	Mouthpiece	5+ Years	

Patent Number : CN 219049892U



LB2124

Chamber Volume: 155ml
Square Mouthpiece

Cat. NO.	Size	Applicable People	Color
LB212401	S	<18 Months	Orange, Yellow, Dark blue, Green, Whithe
LB212402	M	1-5 Years	
LB212403	L	5+ Years	
LB212404	Mouthpiece	5+ Years	

Face Mask (Spacer)

The masks is fit for universal inhaler



Silicone Mask Square	Silicone Mask Round	PVC Mask Square	PVC Mask Round	Size	Application People
LB-SS-S	LB-SR-S	LB-PS-S	LB-PR-S	S	0-18 Months
LB-SS-M	LB-SR-M	LB-PS-M	LB-PR-M	M	1-5 Years
LB-SS-L	LB-SR-L	LB-PS-L	LB-PR-L	L	5+ Years



LB1810 Three balls Incentive Spirometer

The three balls Incentive Spirometer is designed to encourage the user to take deep breath. The spirometer provides patient feedback by showing how well the patient is expanding their lungs. This deep breathing encourages the clearance of phlegm from the lungs.

Cat. NO.	Size
LB181001	600cc/sec+900cc/sec+1200cc/sec
LB181002	700cc/sec+1300cc/sec+1900cc/sec




LB1811 Volumetric Incentive Spirometer

Size: 3000ml, 5000ml

Incentive Spirometer is handheld medical device that measures the volume of patients' breath. It helps patients' lungs recover after surgery or lung illness, keeping them active and free of fluid. The patients are encouraged to breath in from the device as slowly and as deeply as possible.

Cat. NO.	Size
LB181101	3000cc
LB181102	5000cc

Patent Number : CN 213789846U



LB1111 Pediatric Cartoon Oxygen Mask

Made of clear, medical grade PVC.
Elastic strap with/without latex.
Medium concentration, approx 50%.
Length- with 2m tubing, different length can be customized.
Anti-crush tube can ensure oxygen flow even if the tube is kinked.
Packing: packed in individual PE pack, 50 pcs/carton or 100 pcs/carton.

Code No.	Pattern	Pattern Color
LB1111P	Panda	Black, Blue, Green etc for optional
LB1111E	Elephant	Orange, Blue, Green etc for optional



LB2111 Pediatric Cartoon Nebulizer Mask with 360° swivel connector

Without aluminum nose clip can fit for MRI
Made of clear, medical grade PVC
Swivel 360° connector
Without adjustable nose clip, elastic strap with/without latex
Length: with 2m tubing, different length can be customized,
Anti-crush tube can ensure oxygen follow even if the tube is kinked.
Packing: packed in individual PE pack. 50 pcs/carton or 100 pcs/carton.

Code No.	Pattern	Pattern Color
LB2111P	Panda	Black, Blue, Green etc for optional
LB2111E	Elephant	Orange, Blue, Green etc for optional



LB1110 Simple Oxygen Mask

Made of clear, medical grade PVC
Adjustable nose clip
Anti-crush tube

Cat. NO.	Size
LB111000	XS, Neonate
LB111001	S, Pediatric Standard
LB111002	M, Pediatric Elongated
LB111003	L, Adult Asian
LB111004	L+, Adult European
LB111005	XL, Adult Elongated



LB1210 Venturi Mask

Made of clear, medical grade PVC
Two color coded diluters for seven concentrations
Green diluter for 24%, 26%, 28%, 30%
White diluter for 35%, 40%, 50%
Adjustable nose clip
Anti-crush tube

Cat. NO.	Size
LB121000	XS, Neonate
LB121001	S, Pediatric Standard
LB121002	M, Pediatric Elongated
LB121003	L, Adult Asian
LB121004	L+, Adult European
LB121005	XL, Adult Elongated



LB1220 Adjustable Venturi Mask

Made of clear, medical grade PVC
Single diluter for one concentration, seven diluters can be selected
Adjustable nose clip, Adaptor for high humidity, Anti-crush tube
With seven diluters 24% (blue) 28% (yellow) 31% (white) 35% (green) 40% (pink) 50% (orange) 60% (red)


Cat. NO.		Size
6 diluters	7 diluters	
LB122000A	LB122000B	XS, Neonate
LB122001A	LB122001B	S, Pediatric Standard
LB122002A	LB122002B	M, Pediatric Elongated
LB122003A	LB122003B	L, Adult Asian
LB122004A	LB122004B	L+, Adult European
LB122005A	LB122005B	XL, Adult Elongated



LB1310 Oxygen Mask with Reservoir Bag

Made of clear, medical grade PVC
Standard connector, non-rebreathing
With 1000ml/600ml reservoir bag
Anti-crush tube

Cat. NO.	Size
LB131000	XS, Neonate
LB131001	S, Pediatric Standard
LB131002	M, Pediatric Elongated
LB131003	L, Adult Asian
LB131004	L+, Adult European
LB131005	XL, Adult Elongated





Nasal Oxygen Cannula

Made of clear, medical grade PVC
Dipped prong, straight prong, curved prong is available
Green transparent, white transparent and light blue transparent colors are available.
Different lengths can be customized.
Anti-crush tube can ensure oxygen flow even if the tube is kinked.
Ear Guards type is available.

Size	Cat.No.: LB1411 (With Curved Prong)	Cat.No.: LB1412 (With Straight Prong)	Cat.No.: LB1413 (With Dipped Prong)
XXS (Neonate)	/	LB141201	LB141301
XS (Infant)	/	LB141202	LB141302
S (Child)	LB141103	LB141203	LB141303
L (Adult)	LB141104	LB141204	LB141304

**LB1610
Oxygen Connecting Tube**

Made of clear, medical grade PVC
The length of the tube can be customized.



**LB1510
Tracheostomy Mask**

Made of clear, medical grade PVC
With elastic strap with/without latex
Swivel 360° adaptor with 18F/22M port.
Excellent compatible with corrugated tubing & other connection device.

Cat. NO.	Description
LB151001	Pediatric Mask
LB151002	Adult Mask
LB151101	Pediatric Mask, Tubing Connector, 2.0m Tubing
LB151102	Adult Mask, Tubing Connector, 2.0m Tubing
LB151201	Pediatric Mask, Button Type
LB151202	Adult Mask, Button Type
LB151301	Pediatric Mask, Button Type, Tubing Connector, 2.0m Tubing
LB151302	Adult Mask, Button Type, Tubing Connector, 2.0m Tubing

Nebulizer Jar



**LB20A
Fixed Nebulizer Jar**

Volume: 6cc
Under the condition of gas flow 7L/min, atomizing volume about 200mg/min, the percentage of particulate with the diameter less than 5um is more than 70% in every 200mg.



**LB20B
Fixed Nebulizer Jar**

Volume: 6cc
Under the condition of gas flow 7L/min, atomizing volume about 280mg/min, the percentage of particulate with the diameter less than 5um is more than 70% in every 280mg.

Patent Number :
ZL 2014 2 0766313.3



**LB20C
Fixed Nebulizer Jar**

Volume: 10cc
Under the condition of gas flow 7L/min, atomizing volume about 250mg/min, the percentage of particulate with the diameter less than 5um is more than 70% in every 250mg.

Patent Number :
ZL 201520637654.5



**LB20D
Adjustable Nebulizer Jar**

Volume: 6cc
Under the condition of gas flow 7L/min, atomizing volume about 150-450mg/min, the percentage of particulate with the diameter less than 5um is more than 70%.

Patent Number :
ZL 201520289256.9



**LB20E
Adjustable Nebulizer Jar**

Volume: 10cc
Under the condition of gas flow 7L/min, atomizing volume about 120-450mg/min, the percentage of particulate with the diameter less than 5um is more than 70%.

Patent Number :
CN 209967301U



**LB20F
Adjustable Nebulizer Jar**

Volume: 10cc
Adjustable Nebulizer Jar
Under the condition of gas flow 7L/min, atomizing volume about 120-450mg/min, the percentage of particulate with the diameter less than 5um is more than 71%.

Patent Number :
CN 209967301U



LB2011 Nebulizer Mask

Anti-crush tube can ensure oxygen flow even if the tube is kinked. Different types of nebulizer jars are available.

Cat. NO.	Size
LB201100	XS, Neonate
LB201101	S, Pediatric Standard
LB201102	M, Pediatric Elongated
LB201103	L, Adult Asian
LB201104	L+, Adult European
LB201105	XL, Adult Elongated



LB2021 Nebulizer Mask with 360° swivel connector

It is without aluminum nose clip, and can fit for MRI. Anti-crush tube can ensure oxygen flow even if the tube is kinked. Different types of nebulizer jars are available. Size: Adult, Child



LB2012 Nebulizer Kit

Green transparent, white transparent and light blue transparent colors are available. Anti-crush tube can ensure oxygen flow even if the tube is kinked. Different types of nebulizer jars are available.

Code No.	Description
LB201201	Adult Aerosol Mask, Pediatric Aerosol Mask, Nebulizer Chamber, T-piece, Mouthpiece, Corrugated Tube, 2.0m Tubing
LB201202	Adult Aerosol Mask, Pediatric Aerosol Mask, Nebulizer Chamber, T-piece, Mouthpiece, 2.0m Tubing
LB201203	Adult Aerosol Mask, Pediatric Aerosol Mask, Nebulizer Chamber, Angled Mouthpiece, 2.0m Tubing
LB201204	Adult Aerosol Mask, Pediatric Aerosol Mask, Nebulizer Chamber, 2.0m Tubing
LB2013	Nebulizer Chamber, T-piece, Corrugated tube, Mouthpiece, 2.0m Tubing
LB2014	Nebulizer Chamber, T-piece, Mouthpiece, 2.0m Tubing
LB2015	Nebulizer Chamber, Angled Mouthpiece, 2.0m Tubing




LB2032 Hand Holding Nebulizer Set

Made of clear, medical PVC
 With 1pc adult mask and 1pc child mask mouthpiece,
 2m oxygen tubing
 Volume: 10cc



LB60 Suction Connecting Tubing

- 1) Made of transparent material for better visualization
- 2) Striated walls of tube provide superior strength and anti-kinking
- 3) Multiple choices of length
- 4) Size: 1/4", 3/16", 9/32"
- 5) Tube Length Customized



Suction Yankauer Handle

Size	Type
LB6101	Crown tip with vent
LB6102	Crown tip without vent
LB6103	Flat tip with vent
LB6104	Flat tip without vent



LB61 Suction Connecting Tubing with Yankauer Handle

1. Made of transparent material
2. Can be combined with different types of yankauer handle and connecting tube to suit various requirement
3. Multiple choices of lengths
4. Can choose different connectors



**LB3011
Anesthesia Mask**

Made from Latex/DEHP free PVC material.
Air cushion assures comfortable face fitting.
Tail check valve/Upright check valve without check valve type are available.
Single Patient use.
Color-coded hook options.

Cat.No.	Size	Color	Connector
LB301111	Size 000#	/	15mm OD
LB301110	Size 00#	/	15mm OD
LB301100	Size 0# - Neonate	/	15mm OD
LB301101	Size 1# - Infant	/	15mm OD
LB301102	Size 2# - Pediatric	Yellow	22mm ID
LB301103	Size 3# - Adult-small	Green	22mm ID
LB301104	Size 4# - Adult-medium	Red	22mm ID
LB301105	Size 5# - Adult-large	Blue	22mm ID
LB301106	Size 6# - Adult-Extra large	Orange	22mm ID



**LB3021
Soft Anesthesia Mask**

Made from Latex/DEHP free PVC material.
Soft adjustable face cushion assures comfortable face fitting.
With Upright Check valve.
Single Patient use.
No scents or with Vanilla scents are available.

Cat.No.	Size	Color	Connector
LB302101	Size 0# - Premature	White	15mm OD
LB302102	Size 1# - Neonate	Pink	15mm OD
LB302103	Size 2# - Infant	Yellow	15mm OD
LB302104	Size 3# - Pediatric	Green	22mm ID
LB302105	Size 4# - Adult-S	Red	22mm ID
LB302106	Size 5# - Adult-M	Blue	22mm ID
LB302107	Size 6# - Adult-L	Orange	22mm ID



**LB3012
PVC Free Anesthesia Mask**

Made from TPE and PP Material, contain no PVC and are DEHP FREE
The products provides a comfortable superior airtight fit with a clear shell to enhance patient face visibility.
This products ancillary use with anesthesia breathing circuit.
Non-reusable, single use only

Cat.No.	Size	Connector
LB301201	Size 1# - Neonate	15mm OD
LB301202	Size 2# - Infant	15mm OD
LB301203	Size 3# - Pediatric	22mm ID
LB301204	Size 4# - Adult-S	22mm ID
LB301205	Size 5# - Adult-M	22mm ID
LB301206	Size 6# - Adult-L	22mm ID



**LB3030
Silicone Anesthesia Mask (One-piece)**

All-in-one piece soft silicone cushion for better seal and comfort.
Ergonomic shape for lower dead space and better air circulation.
Made of 100% clear, medical grade silicon, latex free
Reusable, all parts are fully autoclavable up to 134°C

Cat.No.	Size	Connector
LB303000	Size 0# - Neonate	15mm OD
LB303001	Size 1# - Infant	15mm OD
LB303002	Size 2# - Pediatric	22mm ID
LB303003	Size 3# - Adult-S	22mm ID
LB303004	Size 4# - Adult-M	22mm ID
LB303005	Size 5# - Adult-L	22mm ID
LB303006	Size 6# - Adult-XL	22mm ID



LB3031 Silicone Anesthesia Mask (rounded)

All-in-one piece soft silicone cushion for better seal and comfort.
Ergonomic shape for lower dead space and better air circulation.
Made of 100% clear, medical grade silicon, latex free
Reusable, all parts are fully autoclavable up to 134°C

Cat.No.	Size	Connector
LB303130	000# - Premature	15mm OD
LB303120	00# - Premature	15mm OD
LB303110	0# - Neonate	15mm OD
LB303101	1# - Infant	15mm OD
LB303102	2# - Child	22mm ID



**LB3035
Endoscope Mask**

Features
1. Used in painless gastroscopy
2. Intubation and ventilation can be operated and reduce the operation time
3. Extend joint can connect oxygen plant
4. Various specification fiberscope apertures are convenient for work-up and treatment

Specification

Mask Size	Hole Diameter
3#	10mm
4#	12mm
5#	10mm, 12mm



LB3032 Silicone Anesthesia Mask (two-piece)

Soft silicone cushion encased in hard shell for ease to handle and better seal.
Made of 100% clear, medical grade silicon and pc material, latex free.
This products ancillary use with disposable anesthesia breathing circuit.
Reusable, all parts are fully autoclavable up to 134°C

Cat.No.	Size	Connector
LB303203	3# - Adult-Small	22 mm ID
LB303204	4# - Adult-Medium	22 mm ID
LB303205	5# - Adult-Large	22 mm ID



**LB3041
Silicone CPAP Mask and Headgear**

Type: Full mask and nasal mask CPAP mask are available
Size: L,M,S
Simple and comfortable adjustable masks that assist patients equipped with CPAP(Continuous Positive Airway Pressure) or NIV(Non Invasive Ventilation) ventilators.
Headgear quick clips allow for fast and convenient removal of the mask
The ultra-flexible forehead support reduces contact with your skin.

Anesthesia Breathing Circuit

Anesthesia breathing circuit is bridge of anesthesia machine with patient, input anesthesia gases and output expiratory air from the patient, to accomplish the normal exchange of oxygen and carbon dioxide gas.

Disposable anesthesia circuit can prevent patients from infecting, establish manual airway and match with anesthesia machine and respirator. Non-toxic, insipid, made of medical degree material, flexible transparent, air tightness.



LB431 Anesthesia Circuit-Normal (Adult/Pediatric/Infant)



LB431E Anesthesia Circuit-Extendable (Adult/Pediatric)




LB431R Anesthesia Circuit-Reinforced/Smoothbore (Adult/Pediatric)



LB4330 Anesthesia Circuit-Coaxial (Adult)



LB432 Breathing Circuit-Corrugated (Adult/Pediatric)



LB432R Breathing Circuit-Reinforced/Smoothbore (Adult/Pediatric)



Anesthesia Circuit Kit

Cat.No.	Circuit Type	Size
LB4510	Corrugated	Adult, Pediatric, Infant
LB4511	Extendable	Adult, Pediatric
LB4512	Smoothbore	Adult, Pediatric

Optional Items

Optional Items	Available or Regular Specification
Anesthesia Breathing Cuff	Latex-free, Latex, Silicone
2-Way connector	22M-22M
Filter	HMEF, BV
CO ₂ Sampling line	3m
Extra Limb	0.4m
Pressure Monitoring Line	3m
Mask	PVC, Soft PVC, Silicone
Catheter Mount	Normal, Extendable, Reinforced



LB4410R
LB4410E
LB4410



LB482002 Y Type
LB482003 U Type

LB4410 Catheter Mount (Normal/Extendable/Smoothbore) LB4820 Water Trap


Anesthesia Breathing Bag



0.5L 1L 2L 3L



0.5L 1L 2L 3L



0.5L 1L 2L 3L

LB4910

Anesthesia Breathing Bag-Latex Free

Cat.No.	Size
LB491000	0.5L
LB491001	1L
LB491002	2L
LB491003	3L

LB4920

Anesthesia Breathing Bag-Latex

Cat.No.	Size
LB492000	0.5L
LB492001	1L
LB492002	2L
LB492003	3L

LB4930

Anesthesia Reservoir - Silicone

Cat.No.	Size
LB493000	0.5L
LB493001	1L
LB493002	2L
LB493003	3L



LB1440
High Flow Nasal Cannula

Cat.No.	Description	Package
LB144001	High Flow Nasal Cannula, S	50pcs/Box
LB144002	High Flow Nasal Cannula, M	50pcs/Box
LB144003	High Flow Nasal Cannula, L	50pcs/Box



LB4610
Humidifier Chambers

Heating Plate Diameter: 12cm



Type B

Filter

Disposable breathing filter made of medical degree material, mainly used for filtering bacteria and increase gas humidity degree, effectively prevent patient infecting and mainly used as medical accessory. Efficiency filtration, humidity exchange function.




LB411001 BV Filter



LB411002 BV Filter



LB411003 BV Filter



LB421001 HME Filter



LB421002 HME Filter



LB421004 HME Filter



LB421005 HME Filter



LB421006 HME Filter

Resuscitator

Optional Accessories: Guedel Airway, Manometer, Peep Valve, B-type Intake Valve, etc.



LB3110
Silicone Resuscitator

Built-in 40-60cm H₂O pressure relief valve.
The set includes oxygen reservoir bag, silicone mask and oxygen tubing.
100% Medical level silicone material.
Latex-free components.
All-in-one intake valve connects to oxygen reservoir bag directly.
Additional accessories and customized label are available.
Built-in 30MM exhale port for PEEP or filter attachment.
All parts (excluding the oxygen tubing & oxygen reservoir) are fully autoclavable up to 134°C.



LB3120
SEBS Resuscitator

Single patient use only.
Built-in 60cm H₂O pressure relief valve.
The set includes oxygen reservoir bag, PVC and oxygen tubing.
100% Medical level SEBS material.
Latex-free components.
All-in-one intake valve connects to oxygen reservoir bag directly.
Additional accessories and customized label are available.
Built-in 30MM exhale port for PEEP or filter attachment.



LB3110L
Liquid Silicone Resuscitator

Built-in 40-60cm H₂O pressure relief valve.
The set includes oxygen reservoir bag, silicone mask and oxygen tubing.
100% Medical level silicone material.
Latex-free components.
All-in-one intake valve connects to oxygen reservoir bag directly.
Additional accessories and customized label are available.
Built-in 30MM exhale port for PEEP or filter attachment.
All parts (excluding the oxygen tubing & oxygen reservoir) are fully autoclavable up to 134°C.
Transparency and resilience is better.



LB3130
PVC Resuscitator

Single patient use only.
Built-in 60cm H₂O pressure relief valve.
The set includes oxygen reservoir bag, PVC mask and oxygen tubing.
100% Medical level PVC material.
Latex-free components.
All-in-one intake valve connects to oxygen reservoir bag directly.
Additional accessories and customized label are available.
Built-in 30MM exhale port for PEEP or filter attachment.



LB3510 CPR Mask

Features:
 Non-reusable, single use only.
 Foldable and non-foldable mask design.
 High transparency for better visualization.
 100% Medical level PVC material.
 Latex-free components.

Specification:
 Foldable cushion mask with oxygen port.
 One way valve with happy efficiency filter.
 Elastic head trap.
 Individual colorful hard case package



LB3220 Reusable Silicone Laryngeal Mask

Reusable Silicone Laryngeal Mask are designed to be sterilized and used up to 40 times
 Durable high temperature resistant material suitable for autoclave sterilization
 Additional user card for record of device usage times and patient information is supplied in each pouch

Cat. NO. (Standard)	Cat. NO. (Reinforced)	Type	Applicable People
LB322010	LB322010R	1.0	< 5kg
LB322015	LB322015R	1.5	5-10kg
LB322020	LB322020R	2.0	10-20kg
LB322025	LB322025R	2.5	20-30kg
LB322030	LB322030R	3.0	30-50kg
LB322040	LB322040R	4.0	50-70kg
LB322050	LB322050R	5.0	> 70kg



LB3210 Disposable Silicone Laryngeal Mask

Medical grade silicone for super compatibility and patient comfort
 Available with seven sizes for all patient use

Cat. NO. (Standard)	Cat. NO. (Reinforced)	Type	Applicable People
LB321010	LB321010R	1.0	< 5kg
LB321015	LB321015R	1.5	5-10kg
LB321020	LB321020R	2.0	10-20kg
LB321025	LB321025R	2.5	20-30kg
LB321030	LB321030R	3.0	30-50kg
LB321040	LB321040R	4.0	50-70kg
LB321050	LB321050R	5.0	> 70kg



LB3230 Disposable PVC Laryngeal Mask

Medical grade PVC with good compatibility
 Available with seven sizes for all patient use

Cat. NO. (Standard)	Cat. NO. (Reinforced)	Type	Applicable People
LB323010	LB323010R	1.0	< 5kg
LB323015	LB323015R	1.5	5-10kg
LB323020	LB323020R	2.0	10-20kg
LB323025	LB323025R	2.5	20-30kg
LB323030	LB323030R	3.0	30-50kg
LB323040	LB323040R	4.0	50-70kg
LB323050	LB323050R	5.0	> 70kg

Oropharyngeal Airways

Smooth integrated design for optimum patient comfort and safety.
 Half hard non-toxic polyethylene material with color code.

LB3310 Guedel Type

Cat. NO.	Size, Length	Color Code
LB3310040	000#, 40mm	Pink
LB3310050	00#, 50mm	Blue
LB3310060	0#, 60mm	Black
LB3310070	1#, 70mm	White
LB3310080	2#, 80mm	Green
LB3310090	3#, 90mm	Yellow
LB3310100	4#, 100mm	Red
LB3310110	5#, 110mm	Orange
LB3310120	6#, 120mm	Purple



LB3320 Berman Type

Cat. NO.	Size, Length
LB3320040	000#, 40mm
LB3320050	00#, 50 mm
LB3320060	0#, 60 mm
LB3320070	1#, 70 mm
LB3320080	2#, 80 mm
LB3320090	3#, 90 mm
LB3320100	4#, 100mm
LB3320110	5#, 110 mm



LB3410 Nasopharyngeal Airway

Made from non-toxic PVC, soft and comfortable
 Only for nasal airway

Cat. NO.	Size ID(mm)
LB3410030	3.0 mm
LB3410035	3.5 mm
LB3410040	4.0 mm
LB3410045	4.5 mm
LB3410050	5.0 mm
LB3410055	5.5 mm
LB3410060	6.0 mm
LB3410065	6.5 mm
LB3410070	7.0 mm
LB3410075	7.5 mm
LB3410080	8.0 mm
LB3410085	8.5 mm
LB3410090	9.0 mm
LB3410095	9.5 mm



Endotracheal Tube

X-ray opaque line throughout the length of the tube
Non-toxic PVC, transparent, soft and smooth

Standard Endotracheal Tube

LB5010 Uncuffed

Cat. NO.	Size ID (mm)	Cat. NO.	Size ID (mm)
LB5010020	2.0	LB5010065	6.5
LB5010025	2.5	LB5010070	7.0
LB5010030	3.0	LB5010075	7.5
LB5010035	3.5	LB5010080	8.0
LB5010040	4.0	LB5010085	8.5
LB5010045	4.5	LB5010090	9.0
LB5010050	5.0	LB5010095	9.5
LB5010055	5.5	LB5010100	10.0
LB5010060	6.0		



LB5010C High volume low pressure cuff

Cat. NO.	Size ID (mm)	Cat. NO.	Size ID (mm)
LB5010025C	2.5	LB5010065C	6.5
LB5010030C	3.0	LB5010070C	7.0
LB5010035C	3.5	LB5010075C	7.5
LB5010040C	4.0	LB5010080C	8.0
LB5010045C	4.5	LB5010085C	8.5
LB5010050C	5.0	LB5010090C	9.0
LB5010055C	5.5	LB5010095C	9.5
LB5010060C	6.0	LB5010100C	10.0



Oral Preformed Endotracheal Tube

LB5020 Uncuffed

Cat. NO.	Size ID (mm)	Cat. NO.	Size ID (mm)
LB5020020	2.0	LB5020065	6.5
LB5020025	2.5	LB5020070	7.0
LB5020030	3.0	LB5020075	7.5
LB5020035	3.5	LB5020080	8.0
LB5020040	4.0	LB5020085	8.5
LB5020045	4.5	LB5020090	9.0
LB5020050	5.0	LB5020095	9.5
LB5020055	5.5	LB5020100	10.0
LB5020060	6.0		



LB5020C High Volume Low Pressure Cuff

Cat. NO.	Size ID (mm)	Cat. NO.	Size ID (mm)
LB5020030C	3.0	LB5020070C	7.0
LB5020035C	3.5	LB5020075C	7.5
LB5020040C	4.0	LB5020080C	8.0
LB5020045C	4.5	LB5020085C	8.5
LB5020050C	5.0	LB5020090C	9.0
LB5020055C	5.5	LB5020095C	9.5
LB5020060C	6.0	LB5020100C	10.0
LB5020065C	6.5		



Endotracheal Tube

X-ray opaque line throughout the length of the tube
Non-toxic PVC, transparent, soft and smooth

Nasal Preformed Endotracheal Tube

LB5030 Uncuffed

Cat. NO.	Size ID (mm)	Cat. NO.	Size ID (mm)
LB5030020	2.0	LB5030065	6.5
LB5030030	3.0	LB5030070	7.0
LB5030035	3.5	LB5030075	7.5
LB5030040	4.0	LB5030080	8.0
LB5030045	4.5	LB5030085	8.5
LB5030050	5.0	LB5030090	9.0
LB5030055	5.5	LB5030095	9.5
LB5030060	6.0	LB5030100	10.0



LB5030C High Volume Low Pressure Cuff

Cat. NO.	Size ID (mm)	Cat. NO.	Size ID (mm)
LB5030030C	3.0	LB5030065C	6.5
LB5030035C	3.5	LB5030070C	7.0
LB5030040C	4.0	LB5030075C	7.5
LB5030045C	4.5	LB5030080C	8.0
LB5030050C	5.0	LB5030085C	8.5
LB5030055C	5.5	LB5030090C	9.0
LB5030060C	6.0	LB5030095C	9.5
		LB5030100C	10.0



Reinforced Endotracheal Tube

LB5040 Uncuffed

Cat. NO.	Size ID (mm)	Cat. NO.	Size ID (mm)
LB5040030	3.0	LB5040065	6.5
LB5040035	3.5	LB5040070	7.0
LB5040040	4.0	LB5040075	7.5
LB5040045	4.5	LB5040080	8.0
LB5040050	5.0	LB5040085	8.5
LB5040055	5.5	LB5040090	9.0
LB5040060	6.0	LB5040095	9.5



LB5040C High Volume Low Pressure Cuff

Cat. NO.	Size ID (mm)	Cat. NO.	Size ID (mm)
LB5040030C	3.0	LB5040065C	6.5
LB5040035C	3.5	LB5040070C	7.0
LB5040040C	4.0	LB5040075C	7.5
LB5040045C	4.5	LB5040080C	8.0
LB5040050C	5.0	LB5040085C	8.5
LB5040055C	5.5	LB5040090C	9.0
LB5040060C	6.0	LB5040095C	9.5



LB5011C Endotracheal Tube with Evacuation Lumen

Cat. NO.	Size ID (mm)	Cat. NO.	Size ID (mm)
LB5011050C	5.0	LB5011085C	8.5
LB5011055C	5.5	LB5011090C	9.0
LB5011060C	6.0		
LB5011065C	6.5		
LB5011070C	7.0		
LB5011075C	7.5		
LB5011080C	8.0		



LB5050 Double-Lumen Endobronchial Tube

It is used for one-lung ventilation
 Made of PVC in medical grade
 Available with Left-sided and Right-sided types
 The cuffs, two lumens and pilot balloons are imprinted with Tracheal (clear) or Bronchial (blue) for easy identification
 Smooth inner lumens allow easy passage of bronchoscopes
 Pre-loaded with stylet to help tube maintain shape for easier intubation
 Packed with switch connector and two low friction suction catheters

Right-sided		Left-sided	
Cat.No.	Size(Fr)	Cat.No.	Size(Fr)
LB505028R	28	LB505028L	28
LB505032R	32	LB505032L	32
LB505035R	35	LB505035L	35
LB505037R	37	LB505037L	37
LB505039R	39	LB505039L	39
LB505041R	41	LB505041L	41



Endotracheal Tube Introducer (Bougie)

Smooth elbow design of tip can assist intubation decreasing the damage of soft tissue
 scale mark indicates distance to the tip
 recommend for difficult intubation and replacing an endotracheal tube
 Available with solid type and hollow type.
 Hollow type can facilitate delivery of oxygen in emergency


Cat. NO.	Size (Fr)	Type	Applicable ETT
LB521006	6	Solid	< 3.5mm
LB521010	10	Solid	4.0-5.0mm
LB521015	15	Solid	> 5.5mm
LB521010H	10	Hollow	4.0-5.0mm
LB521015H	15	Hollow	> 5.5mm



LB5310 Intubating Stylet

Made of malleable aluminum covered with plastic sheath
 Soft rounded tip for minimizing patient trauma during intubation.
 Flexible use for assistance of intubation.

Cat. NO.	Size (Fr)	Applicable ETT
LB531006	6	2.5-4.5mm
LB531010	10	4.0-6.0mm
LB531012	12	5.5-6.5mm
LB531014	14	6.0-10.0mm



Tracheostomy Tube

LB5110 Uncuffed

Cat. NO.	Size ID (mm)	Cat. NO.	Size ID (mm)
LB5110030	3.0	LB5110070	7.0
LB5110035	3.5	LB5110075	7.5
LB5110040	4.0	LB5110080	8.0
LB5110045	4.5	LB5110085	8.5
LB5110050	5.0	LB5110090	9.0
LB5110055	5.5	LB5110095	9.5
LB5110060	6.0	LB5110100	10.0
LB5110065	6.5		



LB5110C High Volume Low Pressure Cuff

Cat. NO.	Size ID (mm)	Cat. NO.	Size ID (mm)
LB5110030C	3.0	LB5110070C	7.0
LB5110035C	3.5	LB5110075C	7.5
LB5110040C	4.0	LB5110080C	8.0
LB5110045C	4.5	LB5110085C	8.5
LB5110050C	5.0	LB5110090C	9.0
LB5110055C	5.5	LB5110095C	9.5
LB5110060C	6.0	LB5110100C	10.0
LB5110065C	6.5		



LB5110R Reinforced Tracheostomy Tube

Made of medical grade PVC, high volume pressure.
 Rounded and smooth tip of the obturator minimizes tissue injury during intubation.
 Tube with spiral reinforced wire reducing risk of kinking.

Cat. NO.	Size ID (mm)	Cat. NO.	Size ID (mm)
LB5110030R	3.0	LB5110065R	6.5
LB5110035R	3.5	LB5110070R	7.0
LB5110040R	4.0	LB5110075R	7.5
LB5110045R	4.5	LB5110080R	8.0
LB5110050R	5.0	LB5110085R	8.5
LB5110055R	5.5	LB5110090R	9.0
LB5110060R	6.0	LB5110095R	9.5
		LB5110100R	10.0



LB5051 Endobronchial Blocker Tube

Endobronchial blocker tube is a device which can be inserted down a tracheal tube after tracheal intubation so as to block off the right or left main bronchus of the lungs in order to be able to achieve a controlled one lung ventilation. The lung tissue distal to the obstruction will collapse, thus allowing the surgeon's view and access to relevant structures within the thoracic cavity.
 Our blocker tubes are used to achieve lung separation and one lung ventilation as an alternative to Double-lumen tubes(DLT) and serve as a choice for pediatric patients as well.


Blocker Tube Only		Blocker Tube Kit	
Cat. NO.	Size ID (mm)	Cat. NO.	Size ID (mm)
LB505105	5FR	LB505105S	5FR
LB505107	7FR	LB505107S	7FR
LB505109	9FR	LB505109S	9FR





FY0636 Endotracheal Intubation Set


Specification:
Endotracheal tube with stylet 1pc
Suction catheter FR. 14 with finger controller tip 1pc
Suction catheter FR. 16 with finger controller tip 1pc
Gloves 1pair
Guedel airway 1pc
Mouth pad 1pc
Bite block 1pc
Syringe 5cc 1pc
Gauze swab 2pcs
Suction connecting tube 1pc
Face mask 1pc
Cotton balls with silicone oil 2pcs



LB5510 Tracheostomy Tube Holder

Tracheostomy Tube Holders are designed to provide secure positioning and to minimize movement of the tracheostomy tube. Unnecessary movement of the trach tube can result in accidental dislodgement or displacement of the tube, tracheoesophageal fistula, tracheal stenosis, or airway granuloma.
Features:
Made of latex-free, foam-padded cotton for patient comfort
Velcro hooks attaching easily to the neck flange
Wide band tapers at hooks to reduce excess secretion accumulation
Soft materials and wide band provide better skin care

Cat.No.	Size	Type
LB551001	Adult :430x25mm	Two pieces
LB551002	Pediatric :250x20mm	Two pieces
LB551003	Adult :450x35mm	One piece



LB6200 Suction Catheter

Made from non-toxic PVC, non-irritant soft PVC

Product Description :
Tube is frosted or smooth
Open distal end with 2 side holes
With X-ray or without X-ray

Choice of connectors :
LB6210 Plain type
LB6220 Cap-cone
LB6230 Y-type
LB6240 T-type
Specification : Fr5-Fr22
Standard Tube Length : 50cm



LB6310 Feeding Tube

Made from non-toxic PVC, non-irritant soft PVC

Product Description :
Tube is frosted or smooth
Closed distal end with 2 side holes
With X-ray or without X-ray
Scale line
Standard Tube Length: 40/50cm



LB6410 Stomach Tube

Made from non-toxic PVC
non-irritant soft PVC or 100% silicone

Product Description :
Tube is frosted or smooth
Closed distal or open distal end with 4 side holes
With X-ray or without X-ray
Scale line


Specification: Fr5-Fr22
Standard Tube Length: 120cm



LB6420 Ryle's Stomach Tube

Made of medical grade PVC
X-ray opaque line throughout the length
Closed distal end with 4 side holes
Length: 125m
Stainless steel balls are sealed in the tip end;
It is easier to indwell the tube at a certain position.

Size ID(Fr)	Cat.No.	Color
Fr6	LB642006	Light Green
Fr8	LB642008	Blue
Fr10	LB642010	Black
Fr12	LB642012	White
Fr14	LB642014	Green
Fr16	LB642016	Orange
Fr18	LB642018	Red
Fr20	LB642020	Yellow
Fr22	LB642022	Purple
Fr24	LB642024	Dark Blue



LB6510 Rectal Tube

Made from non-toxic PVC , non-irritant soft PVC

Product Description :
Tube is frosted or smooth
Closed distal end with 2 side holes

Specification: Fr20-Fr36
Standard Tube Length: 40cm



LB6600 Nelaton Catheter (PVC)

Made from non-toxic PVC , non-irritant soft PVC


Product Description :
Tube is Play frosted or smooth
Closed distal end with 2 side holes

Specification: Fr6-Fr18
Standard Tube Length: LB6610 Male 40cm
LB6620 Female 20cm

LB6122 Extractor Mucus

Made from non-toxic PVC

Capacity: 25ml
Specification: Fr6-Fr14
Standard Tube Length: 35cm



LB6123 Specimen Container

This product is made up of bottle, cap with silicone tube, flat cap and small piece of seal
Size: 40cc, 80cc
Type of silicone tube: short and long



LB6810 Closed Suction Catheter

Sputum suction pipe sealing film produced by imported resin, strong toughness, feel is exquisite and conducive to health care workers
Sputum suction tube body design scale mark line, facilitate throughout the sputum suction effect
Closed suction phlegm tube adopts full closed type design, can be used to sync with breathing machine, at the time of sputum suction oxygen to sustainable, convenient and safe operation, to avoid cross-infection
Specification: Fr6-Fr18
24hours/72 hours



LB6811 Closed Suction Catheter

Sputum suction pipe sealing film produced by imported resin, strong toughness, feel is exquisite and conducive to health care workers
Sputum suction tube body design scale mark line, facilitate throughout the sputum suction effect
Closed suction phlegm tube adopts full closed type design, can be used to sync with breathing machine, at the time of sputum suction oxygen to sustainable, convenient and safe operation, to avoid cross-infection
Specification: Fr6-Fr18
24hours/72 hours



Foyomed[®]
For your smile



DURABLE
Folding proof tube,
Protective conduit



PRECISE SCALE
Precise scale design,
Convenient for clinical diagnosis



EASY TO USE
Large overturning felling design,
it is convenient and sanitary to use



UROLOGY

Urethral Catheter

- Made from Latex, silicone coated
- Color-coded for size identification
- Sterile by EO
- For single use only
- Hard & soft valves are available
- Available with different volume of balloon



FY0101 2-Way Latex Foley Catheter

Cat. NO.	Size (Fr/Ch)	Balloon Capacity (cc)	Color Code	Length(cm)
FY010106	6	3	Light Red	27
FY010108	8	3-5	Black	27
FY010110	10	3-5	Grey	27
FY010112	12	10	White	40
FY010114	14	10	Green	40
FY010116	16	10	Orange	40
FY010118	18	30	Red	40
FY010120	20	30	Yellow	40
FY010122	22	30	Violet	40
FY010124	24	30	Blue	40
FY010126	26	30	Pink	40

FY0103 2-Way Female Foley Catheter



Cat. NO.	Size (Fr/Ch)	Balloon Capacity (cc)	Color Code	Length(cm)
FY010312	12	10	White	26
FY010314	14	10	Green	26
FY010316	16	10	Orange	26
FY010318	18	30	Red	26
FY010320	20	30	Yellow	26
FY010322	22	30	Violet	26

FY0104 2-Way Tiemann Foley Catheter

Cat. NO.	Size (Fr/Ch)	Balloon Capacity (cc)	Color Code	Length(cm)
FY010412	12	10	White	40
FY010414	14	10	Green	40
FY010416	16	10	Orange	40
FY010418	18	10	Red	40
FY010420	20	10	Yellow	40
FY010422	22	10	Violet	40

FY0102 3-Way Latex Foley Catheter

Cat. NO.	Size (Fr/Ch)	Balloon Capacity (cc)	Color Code	Length(cm)
FY010216	16	30	Orange	40
FY010218	18	30	Red	40
FY010220	20	30	Yellow	40
FY010222	22	30	Violet	40
FY010224	24	30	Blue	40
FY010226	26	30	Pink	40



FY0105 2-Way All Silicone Foley Catheter

Cat. NO.	Size (Fr/Ch)	Balloon Capacity (cc)	Color Code	Length(cm)
FY010506	6	1-5	Light Red	31
FY010508	8	3-5	Black	31
FY010510	10	3-5	Grey	31
FY010512	12	10	White	40
FY010514	14	10	Green	40
FY010516	16	10	Orange	40
FY010518	18	30	Red	40
FY010520	20	30	Yellow	40
FY010522	22	30	Violet	40
FY010524	24	30	Blue	40
FY010526	26	30	Pink	40

FY0106 3-Way Silicone Foley Catheter

Cat. NO.	Size (Fr/Ch)	Balloon Capacity (cc)	Color Code	Length(cm)
FY010616	16	30	Orange	40
FY010618	18	30	Red	40
FY010620	20	30	Yellow	40
FY010622	22	30	Violet	40
FY010624	24	30	Blue	40
FY010626	26	30	Pink	40

FY0107 Urethral Catheter (Red Latex)

Cat. NO.	Size (Fr/Ch)	Length(cm)
FY010706	6	27
FY010708	8	27
FY010710	10	27
FY010712	12	40
FY010714	14	40
FY010716	16	40
FY010718	18	40
FY010720	20	40
FY010722	22	40
FY010724	24	40

FY0108 Malecot Catheter

Cat. NO.	Size (Fr/Ch)	Length(cm)
FY010812	12	40
FY010814	14	40
FY010816	16	40
FY010818	18	40
FY010820	20	40
FY010822	22	40
FY010824	24	40
FY010826	26	40
FY010828	28	40
FY010830	30	40



FY0109 Latex Penrose Drainage Tube

Made from natural latex
 Size: 1/4", 5/16", 3/8", 1/2", 5/8", 3/4", 7/8", 1"
 Length: 450mm
 Packed in individual blister packing



FY0109S Silicone Penrose Drainage Tube

100% medical grade silicone for superior biocompatibility
 X-Ray opaque line throughout the length
 Soft and smooth tube minimizes trauma, tooth-like to minimize clog or obstruct
 Clear connecting tube provides inosection and drainage visibility
 Sterilized double package
 Size: 6mm x 30cm, 6mm x 50cm, 8mm x 30cm, 8mm x 50cm, 10mm x 30cm, 10mm x 50cm, 12mm x 30cm, 12mm x 50cm



FY0110 Latex T Tube

Made of natural latex , and silicone coated
 Used for open gallbladder
 The stem is often connected to a bile bag
 Size : Fr10-Fr26



FY0111 Latex Condom Catheter

Made of natural latex
 By connecting pipe, fixture, membrane, membrane set of edge and sticky
 Size : 25mm , 30mm , 35mm , 40mm



FY0121 Silicone Condom Catheter

Self-adhesive.
 Made of 100% imported medical grade silicone
 Pop On: suitable for old person or sports person.
 Standard: suitable for regular length and adhesive for daily wear.
 Wide Band: More adhesive and viscous force than standard for added security.

Cat.No.	Type	Diameter				
		S	M	L	XL	XXL
FY0121P	Pop On	24mm	28mm	31mm	35mm	40mm
FY0121S	Standard	24mm	28mm	31mm	35mm	40mm
FY0121W	Wide Band	24mm	28mm	31mm	35mm	40mm



FY0123 Ureteral Stent

A double J ureteral stent, or ureteric stent, is inserted into the ureter to prevent or treat obstruction of the urine flow from the kidney into the bladder. J shaped curls are present at both ends to hold the tube in place and prevent migration .
 Features:
 Made of medical grade polyurethane (PU) with superior biocompatibility and less encrustation
 Pigtail shapes preventing stent migration
 Smooth surface for easy placement and withdrawal
 Thin wall for maximum drainage
 Clear markings facilitating confirmation of placement
 Good radiopacity for easy visualization
 Available with different designs, coatings, length and diameters

Size : 4Fr, 4.6 Fr, 5 Fr, 6 Fr, 7 Fr, 8 Fr
 Length (cm) : 24-30cm

FY0124 Guide Wire


Guide Wires are used to ensure rapid access, easy advancement and replacement of medical devices during urological procedures .
 Material:Stainless Steel,PercJ,PTFE Guide Wire with Fixed Stainless Steel Core,PTFE coated shaft, Nitinol core, hydrophilic Pebax Tip, PTFE sheath, Nitinol core, black and yellow strips
 Size: 0.028",0.032",0.035",0.038"
 Length (cm): 90cm,150cm



Economic Urine Bag

Capacity: 2000ml
 With/without non-return Valve
 Inlet tube: O.D 6.5mm; 90cm Length
 Sterile or non sterile
 Packing: Individual PE Bag, 10pcs/Inner Bag, 250pcs/Ctn

Cat. NO.	Description
FY0201	Without bottom outlet
FY0202	Push Pull Valve
FY0203	Twist Turn Valve
FY0204	T-Tap Valve



FY0206
Urinary Drainage Bag

Capacity: 2000ml
T-Tap Valve
With ribbon
Anti-Reflux Flutter Valve
Inlet tube: O.D 9mm; 100cm Length
Needleless Sampling Port with Cap
Sterilization: EO
Packing: Individual PE Bag, 100pcs/Ctn



FY0207
Urinary Drainage Bag

Capacity: 2000ml
T-Tap Valve
Anti-Reflux Flutter Valve
Inlet tube: O.D 9mm; 100cm Length
Needleless Sampling Port with Cap
With bed sheet clamp
With universal double-hook
Sterilization: EO
Packing: Individual PE Bag, 100pcs/Ctn



FY0209
Urinary Drainage Bag

Capacity: 2000ml
Twist Turn Valve
With ribbon
Anti-Reflux Flutter Valve
Inlet tube: O.D 9mm; 100cm Length
Adapter with cap
Sterilization: EO
Packing: Individual PE Bag, 100pcs/Ctn



FY0210
Urinary Drainage Bag

Capacity: 2000ml
with bed sheet clamp
Anti-Reflux Flutter Valve, Air-Vent on the Bag
Inlet tube: O.D 9mm; 100cm Length
Needleless Sampling Port with Cap
With Single Hanger and Rope Hanger
Silicone Outlet Tube with Folding Clamp
Sterilization: EO
Packing: Individual Blister Packing, 40pcs/Ctn



FY0215
Urinary Drainage Bag

Capacity: 4000ml/5000ml
T-Tap Outlet Valve, Rigid Dog House
Inlet tube: O.D 10mm; 120cm Length
Needleless Sampling Port with Cap
with bed sheet clamp, with Slide Clamp
Two Parts Anti-Reflux Drip Chamber
air-Vent on the Bag
Universal Double-Hook and Rope Hanger
Sterilization: EO
Packing: Individual Blister Packing, 40pcs/Ctn



FY0218
Urinary Drainage Urine Meter

400cc Two-Chamber Urine meter, 2500ml Urine Bag
Inlet tube: O.D 10mm; 120cm Length
Needleless Sampling Port with Cap
with bed sheet clamp, with Slide Clamp
Universal Double-Hook
T-Tap Outlet Valve
Sterilization: EO
Packing: Individual Blister Packing, 10pcs/Ctn


FY0219
Urinary Drainage Urine Meter

200cc Urinometer with Drip Chamber
2000ml Urine Bag
Inlet tube: O.D 10mm; 120cm Length
Needleless Sampling Port with Cap
with bed sheet clamp, with Slide Clamp
Universal Double-Hook and Rope Hanger
T-Tap Outlet Valve
Sterilization: EO
Packing: Individual Blister Packing, 10pcs/Ctn



FY0220
Urinary Drainage Urine Meter

500cc Urine meter with 2000ml Urine Bag
Inlet tube: O.D 10mm; 120cm Length
Drip Chamber with Air-Vent and Anti-Reflux Valve
Needleless Sampling Port with Cap, Plastic Hanger
bed Sheet Clamp, with Slide Clamp
Large Overtuning Valve
Sterilization: EO
Packing: Individual Blister Packing, 10pcs/Ctn



FY0221
Urinary Drainage Leg Bag

500ml, 750ml
Anti-Reflux Flutter Valve
Twist Valve
Two Pre-Assembled Straps
Sterilization: EO
Packing: Individual PE Bag, 48pcs/Ctn



FY0223
Urinary Drainage Leg Bag

750 ml, non-woven fabric on the back, three cavity
Small and white T valve plug,
Inlet tube: O.D 9mm; 30cm Length
Needless sampling port with cap
Silicone tube 4.5cm
Two straps are available
Sterilization: EO
Packing: Individual Blister Packing, 100 pcs/Ctn



FY0224
Urinary Drainage Leg Bag

750 ml, non-woven fabric on the back
Small and white T valve plug,
Inlet tube: O.D 9mm; 30cm Length
Needless sampling port with cap
Two Pre-assembled Straps
Sterilization: EO
Packing: Individual Blister Packing, 100 pcs/Ctn

FY0229
Pediatric Urinary Collection Bag

Safety Adhesive, PVC Bag, Sponge
With foam or without foam is available
Capacity: 100ml, 200ml
Sterilization: EO
Packing: Individual PE Bag, 100 pcs/Inner Box, 5000Pcs/ctn



FY0235
Closed Wound Drainage System Kit

Made of imported medical grade silicone
This product includes 100cc/200cc Silicone Reservoir
7mm/10mm Flat drains or other round drains
and stainless steel trocar
Can use separately or using Trocar for insertion and
combining with silicone Reservoir for fluids syction
and collection
Used for various surgeries
Latex-free
Sterilization: EO
Packing: Individual Blister Packing, 20pcs/Ctn



FY0230
Closed Wound Drainage System (Spring)

Our closed wound drainage system,including 400cc, 3-spring evacuator,
PVC tubing and "Y" connector, and drain with trocar. The 3-spring
evacuator design is Minimize Collapse-Resistance, Packed in double
liner Tyvek, to make sure mdeical cleanness on O.R.operation.
Capacity: 200ml, 400ml
Silicone or PVC tube, 120cm Length
Sterilization: EO
Packing: Individual Blister Packing, 20pcs/Ctn



FY0231
Closed Wound Drainage System (Hollow)

Our closed wound drainage system,including 400cc, 3-spring evacuator,
PVC tubing and "Y" connector, and drain with trocar. The 3-spring
evacuator design is Minimize Collapse-Resistance, Packed in double
liner Tyvek, to make sure mdeical cleanness on O.R.operation.
Capacity: 200ml, 400ml, 600ml, 800ml
Silicone or PVC tube, 120cm Length
Sterilization: EO
Packing: Individual Blister Packing, 20pcs/Ctn
Size: Fr 7, 10, 12, 14, 15, 16, 18, 19



FY0240
Vacuum Drainage Bottle

The Vacuum Drainage Bottles are intended to aspirate, remove,
or sample body fluids. They are shatterproof, light weight, clear plastic
bottles that provide containment of bio-hazardous fluids.
Vacuum drainage bottle system provides a suction source that help to
drive wound drainage. It consists of a vacuum drainage bottle,
a connecting tube and a drain with trocar.
Size: 400ml, 600ml, 1000ml





FY0241 Disposable Drainage Bottle Equipment

Single bottle, 1600ml
 Double bottle, 2000ml
 Three bottles, 2500ml
 Components:
 By attachment, drive pipe, bottle cap, bottle body, dip pipe, support and suspension hook constitution.



FY0241W Disposable Water Seal Chest Drainage System

It is used in pneumothorax or fluid in thorax caused by traumatic hemopneumothorax and spontaneous pneumothorax , etc.
 Size: 2000mL
 Three Types of Connector: Straight, Reducer, Tee can be chosen.



FY0241D Disposable Dry Seal Chest Drainage System

The dry seal valve technique is designed for added patient security which is perfect for the drainage of fluid and gas in the first aid.
 Have hook and strap, can hang the drainage device next to the bed or carry out.
 Size: 500mL
 Three Types of Connector: Straight, Reducer, Tee can be chosen.

FY0334 EVA TPN (Total Parenteral Nutrition) Bag

- Features:**
- Non-PVC material, DEHP free, can be used to fat emulsion and other fat soluble liquid and medicine.
 - No adverse effects on infants and young children and pregnant women.
 - Flexible tube and bag to avoid kinking and break.

Item	Volume	Details
[Image]	500ml	EVA bag body EVA inlet tube
	1000ml	
	2000ml	
	3000ml	
[Image]	500ml	EVA bag body
	1000ml	
	2000ml	
	3000ml	



FY0234 Enteral Feeding Bag Set

Volume: 1000ml, 1200ml
 Gravity Type/Pump Type
 Enteral Feeding Bag has rigid neck for easy filling and handling.
 Strong,dependable hanging ring.
 Easy-to-read graduations and easy-view translucent bag.
 Bottom exit port allows complete drainage.

Ostomy Care

FY0301 One-Piece Colostomy items

One-piece Close Pouch: need to change everyday.
 One-piece Drainable Pouch: can use 2-3day.



FY030101 One Piece Close Pouch

FY030111 One Piece Drainable Pouch (Twist-tie)



FY030112 One Piece Drainable Pouch (Clamp)

FY030113 One Piece Clostomy Bag (Velcro)

FY030310 One Piece Urostomy Bag


FY0302 Two-Piece Colostomy Items

Two-Piece Close Pouch: Basis Plate can use for several days, Bag need to change everyday.
 Two-Piece Drainable Pouch: both can for several days and basis need to be matched together.



FY030211 Two-piece Drainable Pouch (Clamp) FY030201 Two Pieces Colostomy Bag (Closed)

Two-Piece Barrier



FY030401 Hydrocolloid Adhesive(Embedded)

FY030402 Hydrocolloid Adhesive with Nonwoven Border (Embedded)

FY030403 Hydrocolloid Adhesive (Adhesive)

Costomy bag is available for transparent color.

Suction Liner



FY0247A Wide Mouth Type

Size: 2000ML, 3500ML

Product Characteristic:

- Simplicity
- Safety
- Economical
- With filter or without filter

Indication for Use:

This product is used for collecting the biohazardous fluids which produced in ICU and operation room.

Product Features:

1. Double safety design to ensure no liquid is suctioned back to the tubing and suction unit, effectively prevent cross infection and protect the medical staff safe.
2. It's made of environment-friendly polyethylene. It's safe and non-toxic.
3. Easy to use, can be a single use, also can use with the stander.
4. Available in some specifications: with or without shutoff filter, with or without solidifier, 1000mlml, 1500ml and 2000ml.

Product Advantages:

1. Soft and transparent (Convenient for observation liquid waste)
2. Double safety design (To prevent fluid into the suction systems by mistake)
3. It can be connected with relevant suction device or vacuum device, soft bag is seal and reliable inside and can be withstand the pressure of 60kg when fully filled with liquid (Prevent infections liquid overflow caused by medical personnel).



FY0247B Spring Design Type


Size: 2000ML

FY0247C Classic Type

Size: 1000ML, 1500ML, 2000ML, 3000ML

FY0247E Floating Stop Type

Size: 2000ML, 3500ML



Canister

FY0248 Reusable Canister

Type: FY0248A FY0248B/C FY0248E

Specification: 1000ML 1500ML 2000ML 3000ML 3500ML

Wide mouth type or Straight type for different suction liner. Unibody design with anti-impact material, the canister is high intensity, good transparency and easy to check the volume scale. The T junction air way design is easy for the suction connection.

FY0249 Suction Liner System

FY024901 Suction Liner System 1 Canisters tandem


FY024902 Suction Liner System 2 Canisters tandem

FY024903 Suction Liner System 4 Canisters tandem

Our suction liner system mainly composed of a trolley, canister, a negative pressure gauge, suction liner, connecting tube and adapter. Can be used with single, double, triple, or quadruple bottles. Different negative adapter can be chosen, Such as: American standard, British Standard, German Standard, etc



Suction Liner System



Foyomed[®]
For your smile

HYPODERMIC



FY0601
Disposable Hypodermic Syringe

1ml, 2/3ml, 5ml, 10ml, 20ml, 30ml, 50/60ml
Three-part syringes/Two part syringes (luer slip, special design to prevent plunger slipping out)
Luer slip/luer lock
Latex/latex free piston
W/needle, W/O needle
Sterilization: EO
PE/Blister Individual Package



FY0602
Auto Self-destruction Syringe

0.05ml, 0.1ml, 0.5ml, 1ml
Fixed needle 27G-30G
Sterilization: EO
Blister Individual Package, 200pcs/Inner Box



FY0603
Tubercle Bacillus Syringe

1ml
Luer slip
Latex/latex free piston
PE/Blister Individual Package
Sterilization: EO



FY0604
Insulin Syringe

0.3ml, 0.5ml, 1ml
Fixed needle 27G-31G
PE/Blister Individual Package
Sterilization: EO



FY0604S
Insulin Safety Syringe

Product Specification:
Sterilized, non-toxic, non-pyrogenic, single use only
Blister packing
Sizes available: 0.3ml, 0.5ml, 1ml
Sterilized by Ethylene oxide gas

Benefits:
Safety Design: Protect the healthcare operators by pull the safety cartridge to withdraw the needle into cartridge after injection
Low Dead Space: Fixed needle design with very low residual after injection
Sharp Needle: Tri-beveled stainless cannula tip
Clear Graduation: Accuracy in small dose injection

FY0617
Disposable Destruction Syringe

1ml, 2/3ml, 5ml, 10ml
After injection, when you pull the plunger, the plunger will be broken
Sterilization: EO



FY061801
Safety Clip Syringe

1ml, 2/3ml, 5ml, 10ml, 20ml
Safety syringe with protective cap
Luer lock
After injection, push the clip over the needle and lock it onto the clip base
Sterilization: EO



FY061802
Safety Syringe with Retractable Needle

Size: 3ml, 5ml, 10ml



FY0401
Irrigation Syringe

50ml, 60ml, 80ml
Optional: Adapter Cover
Type: Flat, Bulb, Ring
PE/Blister Individual Package
Sterilization: EO




FY0641
Oral Syringe

Oral syringes are applied to dispense liquid solutions and suspensions and can be used to give almost any medication available as a capsule or tablet as an oral liquid. They are also useful for gradually increasing or decreasing the dose of the medication.

Oral syringes provide accurate and precise measurement of medications allowing for flexibility of dosing to meet your individual needs. Liquid medications are easier to swallow and can be helpful in small children, pets, and people who have trouble swallowing capsules or pills.

Size: 1ml, 3ml, 5ml, 10ml
PE/Blister Individual Package with tip cap or adapter is available





**FY0501
Disposable Infusion Set**

For infusion bags and bottles
20/60drops=1±0.1ml
Plastic piercing device with protective cap, drip chamber, soft tube, flow regulator, air-inlet
Latex or latex-free flash bulb
With or without fluid filter
With or Without needle
Different types of components and parts can be changed as customer's request
For gravity use
EO sterile
PE/Blister individual package



**FY0502
Disposable Infusion Set**

For infusion bags and bottles
20/60drops=1±0.1ml
Plastic piercing device with protective cap, drip chamber, soft tube, flow regulator, air-inlet
With irrigation Y set
With or without fluid filter
With or Without needle
Different types of components and parts can be changed as customer's request
For gravity use
EO sterile
PE/Blister individual package




**FY0503
Disposable Feeding Infusion Set**

20drops=1±0.1ml
Spike without air Vent
50mm drip chamber without filter
Big connector with cap
EO sterile
PE/Blister individual package



**FY0504
Burette Infusion Set**

Burette Capacity: 100ml, 150ml
For infusion bags and bottles
20/60 drops = 1±0.1ml
with litter
with Y Site
Different types of components and parts can be changed as customer's requirements
For gravity use
EO gas sterile
PE/Blister Individual Package



**FY0505
Disposable Blood Transfusion Set**


For blood bags
20 drops = 1±0.1ml
Different types of components and parts can be changed as customer's requirements
For gravity use
EO gas sterile
PE/Blister Individual Package
Different accessories are available



Needless Y Set



Injection Bulb



Y Set



**FY0503F
Disposable Infusion Set with Precision Filter**

Precision solution filter: The filter medium has 5µm, 3µm, 2µm, filtration rate more than 90% for micro material above 5µm, 3µm, 2µm, effectively improve intravenous infusion safety, reduce near-term and long term damage caused by particles to human body.
Safety regulator: Adjusting range is 30mm and from zero to maximum, regulator is accurate and thick, undeformable.
"Types: PVC infusion set, PVC infusion set without DEHP, TPE infusion set, such infusion sets with better safety in turn and can meet different needs."



**FY0502M
Infusion Set with Micro Regulator**

Precision flow regulator for fine-tuning the flow rate in the range of 250ml/h for precise flow control.
The basic type liquid filter has a filtration rate of more than 80% for the particles of 15µm; the filter medium of the precision type liquid filter has a design of 5µm, 3µm, and 2µm, which is the filtration rate of 5µm, 3µm and 2µm in the liquid solution. The filtration rate is not less than 90%. It can effectively improve the safety of intravenous infusion and reduce the near-term damage and long-term damage which caused by particles to the human body.



**FY0125
Pre-filled Saline Flush Syringe**

Different size syringe barrel diameter is exact same as standard 10mL syringe barrel diameter to prevent vein or catheter damage.
Pre-Filled Syringe can be sterilized in two ways:
1. FY0125N, Normal Saline Syringes: Sterile fluid path, can not use in a sterile field
2. FY0125S, Sterile field saline syringes : Externally sterile, can use in a sterile field.
0.9% sodium chloride or glycerin
Size: 3ml, 5ml, 10ml
Sterilization Method: Steam Sterilization



**FY050601
Dental Needle**



FY050602 Dental Needle

Cat. NO.	Size
FY050601	27g x 38mm, 30g x 38mm, etc.
FY050602	18G-25G



**FY0605
Disposable Hypodermic Needle**

Size: 18-30G
Length of needle: 1/2", 5/8", 3/4", 1", etc.
Blister Individual Package
Sterilization: EO

**FY0605S
Sterile Safety Needles For Single Use**

Size: 18G-30G
Needle bevel oriented to the safety cover to facilitate low-angle injections
Easy to use, allow one-handed activation via two ways (thumb or surface)
Audible and tactile click confirms safety activation



**FY0606
Scalp Vein Set**

Size: 18-30G
Length of needle: 1/2", 5/8", 3/4", 1", etc.
Tube length: 30cm
Luer slip/luer lock
PE/Blister Individual Package
Sterilization: EO



**FY0606S
Safety Scalp Vein Set**

Instruction:
It has a system of retraction, which can lock the needle after use. Therefore, it can protect health care professionals from needle stick injuries and infection, as well as preventing the reuse of the product.
Size: 18G, 19G, 20G, 21G, 22G, 23G, 25G, 27G
Needle Length: 3/4"



**FY0616
Insulin Pen Needle**

Ultra fine cannula for maximum patient comfort, patient pain with triple beveling for easy and pain decreases
The needle is fit for Humapen (EliLily, USA)
NovopenIII (Novo Nordisk, Denmark)
Sterile, non-pyrogenic
Sterilization: EO



**FY0701
Disposable Blood Bag**

Type : Single, Double, Triple, Quadruple
Size: 250ml, 350ml, 400ml, 500ml
Needle: 16G or 17G
Anticoagulant (CPDA or CPD SAGM) can be chosen
W or W/O red cell filtration or whole blood filtration is available.
Pre-donation sampling device, needle protective device and other accessories can be customized.
Steam sterilization



**FY0508
Heparin Cap**

Material: PC
Sterilization: EO
Individual Hard Blister Packing, 100pcs/box, 5000pcs/ctn



**FY0510
Disposable Combi Stoppers**

Material: ABS
Sterilization: EO
Color: Red, White, Blue, etc.
Individual Hard Blister Packing, 100pcs/box



**FY0512 Disinfection Luer Cap
(Keep Contamination Away)**

Product Features:
Disinfect completely
The inside alcohol can wet the whole surface of the connectors, kill the 99.99% bacteria in 3 mins.
Automatic swab
The inside sponge with alcohol will swab and clean the connector when screw the cap on the connector, there is no need to do extra cleanness.
Comprehensive Protection
The hard housing and integrity structure design will prevent from damage the connector in order to avoid outer contamination.
International Standard
The standard luer lock design, can be adapted with all the luer male connectors as IV catheter connector, infusion connectors, CVC catheter connector, and other luer connectors.



FY062001

FY0620 Needle Free Connector

Features:

1. Transparent housing material: Polycarbonate or Copolyester.
2. Metal-free & MRI-compatible.
3. Latex-free.
4. Multiple insertion at least 100 times/day.
5. Priming Volume: 0.09ml.
6. Perfect flowing rate: 350 ml/min under one meter water pressure, tested by Baihe Technical Team.

Cat. NO.	Description
FY062001	Needle Free Connector
FY062002	Y-port Valves
FY062003	Tube End Valves
FY062004	Positive Pressure Needle Free Connector



FY062002



FY062003



FY062004



FY0513 Transfer Spike

Material:ABS
Package: Blister Packing/Bulk



FY0511 Extension Tube with Mirco-flow

Flow range: 5ml/h-250ml/h
Regulator Materials: ABS
Tube: PVC
Other accessories: PC & ABS
Extension Tube Length: 45cm is common, others are also available.
Sterilization: EO
Packing Details: Poly-bag per pcs, 200pcs/inner box, 20boxes/carton



FY050702 3-Way Stopcock with Extension Tube

The length of extension tube can be changed as customer's requirements
Different types of components and parts can be changed as customer's requirements
Sterilization: EO
Packing: Individual Blister Packing

Tube Length:
10-25cm, 50pcs/box
30-50cm, 50pcs/box
80cm, 50pcs/box
100-120cm, 40pcs/box

Cat. NO.	Description
FY050701	Three Way Stopcock
FY050702	3-Way Stopcock with Extension Tube
FY050703	Disposable Extension Tube

FY0509 Extension Tubes with Injection Site Adapters

1. This device is designed to meet the needs of general I.V. Therapy, Anaesthesia Cardiovascular,
2. ICU&CCU, Recovery&Oncology.
3. Safe for infusion IV medication including most high risk anti-neoplastics, or aspirating blood without risk of hemolysis.



FY050900
Extension Tube with Normal Needle Free Connector



FY05091
Extension Tube with Positive Pressure Needle Free Connectors



FY05092
Extension tube with Normal Connector




FY05093
T-Extension Tube



FY05094
Screw Extension Tube



FY05095
Mini-Quantity Infusion Extension Tube (OD:2.0mm)



LB0630 Central Venous Catheter

Advanced automatic drilling technology and equipment, a more accurate location and the pore diameter of the pipe, surface smooth and without cutting residues and injury.

Special custom polyurethane fine materials, delicate and smooth tube, welding head end, ultra smooth surface treatment and X line to show the position, make a catheter inserted into the blood vessels in more easy and safe.



FY6002 Hemodialysis Blood Tubing Set

Hemodialysis Blood Tubing Set Consists of arterial and venous lines, used during dialysis and are attached with A.V. fistula needle and dialyzer.

- Features:**
1. Medical grade materials displaying good biocompatibility.
 2. Fluid pathway is smooth inside, no clotting or air bubbles.
 3. Tubes are flexible, kink-resistant.
 4. Color Coding.
 5. Pump segment has excellent anti-fatigue performance.
 6. Luer-lock connectors according to international standard.
 7. Various series, good applicability.
 8. Individual and customized services provided.



FY3007 Peripherally Inserted Central Catheter Kit

1. The connector is round and smooth with PASV valve.
 - A) Valve opens inward for sampling under negative pressure.
 - B) Valve opens outward for infusion under positive pressure.
 - C) Valve remains closed to prevent air embolism, blood reflux and clot under normal pressure
2. The catheter is made from premium silicone, which has excellent biocompatibility.
3. Adjustable length of catheter.
4. Detectable under X ray, the end tip enhances radiopacity.
5. The style of catheter ensures placement more convenient and smooth.
6. Variety of kit configuration, initial guide wire and positive, pressure needless connector is available for choice.



FY0631 CVC Set

- | | |
|-------------------------------|--------------------------|
| Standard Kits include: | Optional Compound |
| 1. Central Venous Catheter | 1. 5ml Syringe |
| 2. Guide-wire | 2. Surgical Gloves |
| 3. Vessel Dilator | 3. Surgical Pledget |
| 4. Clamp | 4. Surgery Sheet |
| 5. Fastener: Catheter Clamp | 5. Surgery Towel |
| 6. Introducer Needle | 6. Sterile Brush |
| 7. Introducer Syringe | 7. Gauze Pad |
| 8. Injection Needle | 8. Suture of Needle |
| 9. Injection Cap | 9. Wound Dressing |
| | 10. Scalpel |



FY0632 Abdomen Drainage Catheter Kit

- | | |
|-------------------------------|--------------------------|
| Standard Kits include: | Optional Compound |
| 1. Drainage Catheter | 1. Drainage Catheter |
| 2. Guide-wire | Standard Kit Accessories |
| 3. Vessel Dilator | 2. 5ml Syringe |
| 4. Introducer Needle | 3. Surgical Gloves |
| 5. Introducer Syringe | 4. Surgical Pledget |
| 6. Injection Needle | 5. Surgery Sheet |
| 7. Injection Cap | 6. Surgery Towel |
| 8. Connecting tube | 7. Sterile Brush |
| | 8. Gauze Pad |
| | 9. Suture of Needle |
| | 10. Wound Dressing |
| | 11. Scalpel |



FY0633 Blood Purification Products

- | | |
|-------------------------------|--------------------------|
| Standard Kits include: | Optional Compound |
| 1. Hemodialysis Catheter | 1. Hemodialysis Catheter |
| 2. Vessel Dilator | Standard Kit Accessories |
| 3. Introducer Needle | 2. 5ml Syringe |
| 4. Syringe | 3. Surgical Gloves |
| 5. Guide-wire | 4. Surgical Pledget |
| 6. Adhesive Wound Dressings | 5. Surgery Sheet |
| 7. Herparin Caps | 6. Surgery Towel |
| 8. Scalpel | 7. Sterile Brush |
| 9. Needle with Suture | 8. Gauze Pad |




FY0635 Disposable Pressure Transducer

Disposable Pressure Transducer is made of medical PC and PVC material.

Operating Parameters:
 Operating Pressure Range: -30mmHg~+300mmHg
 Operating Temperature Range: 15°C ~40°C
 Storage Range: -25°C ~+70°C
 Atmosphere Pressure Range: 56.7kPa~113kPa
 Zero Output: ±25mmHg
 Sensitivity: 5 μ V/V/mmHg ±2%
 Input Impedance: 1200 Ω ~3200 Ω
 Output Impedance: 285 Ω ~315 Ω
 Zero Thermal Drift: ±0.3mmHg/°C
 Zero Drift: ≤2mmHg
 Natural Frequency: 40Hz nominal for a standard kit >200Hz for transducer alone
 Sensitivity Thermal Drift: ±0.1%/°C
 Leak Current: <2 μ A at 230VRMS, 50Hz
 Note: rated voltage 6.00VDC; ambient temperature 25 °C.

- Specification:**
- | | |
|--|---------|
| 130cm pressure tubing with M/F connector | 1piece |
| 30cm pressure tubing with M/F connector | 1piece |
| Disopsable pressure transducer | 1piece |
| Flush device 3 cc/h (30cc/h) | 1piece |
| Male non-vented Cap | 2pieces |
| 150cm single IV set | 1piece |
| 3-way stopcock | 2pieces |



FY0612 A.V. Fistula Needle

It is an irregular connection between an artery and a vein. General Type or Safety Type. Fixed wing or rotating wing is available. Size: 15G, 16G, 17G. EO Sterile.



FY0611 Spinal Needle

Size: 18g~28g x 90mm
 Type: Pencil point and quincke point, can with introducer or not
 Sterilization: EO
 Packing: 50pcs/Box, 1000pcs/ctn



FY0651 Epidural Needle

Size: 15G, 16G, 17G, 18G
 Length: 80/90mm



FY0615 Mini Anesthesia Kit

	Combined Kit FY0615-SE	Epidural Set FY0615-E
Epidural Needle: 16G/18G x 80mm	√	√
Spinal Needle: 25Gx110mm, Pencil Point	√	/
Epidural Catheter: Multiport with mark 0.8/1.0mm≥850mm	√	√
Epidural Filter: 0.22um, hydrophobic membrane, luer lock	√	√
LOR Syringe: 7ml, 10ml available	√	√
Catheter Adaptor: TBA, Secure the catheter completely	√	√

FY0615-SE Combined Kit



FY0615-E Epidural Set

I.V. Cannula

Normal/Safety
 Pen type/With wings/With injection port
 Size: 14G-26G
 PTEF/FEP/PU flex catheter
 Sterilization: EO
 Packing: Hard Blister Packing




FY0607 Pen Type



FY0609 With wings and injection port



FY0608 With Wings



FY0901 Disposable Surgical Blade

Material: Carbon or stainless steel sterile
 Blade shapes #9 to 17 are used in histology units and are all compatible with the associated #3 handles (Blade shape 11 can be also used for gardening with a plastic or metal handle).
 Blade shapes #18 to 36, #60 and #70 are used for major invasive procedures and are all compatible with the associated #4 handles
 Sterilization: Gamma
 Packing: 1pc/Sealed Aluminum Foil Bag, 100pcs/box, 5000pcs/ctn



FY0902 Stitch Cutters

All sizes are available in both carbon or stainless steel
 Designed with a rounded bull nose tip, the Index family of Stitch Cutter Blades are ideal for cutting the wide range of modern suture available today whilst minimising patient discomfort
 The range consists of a short standard stitch cutter with a larger grip. The long handle enables a more aseptic technique allowing good visibility whilst working further from the wound site
 Size: 110mm
 Sterilization: Gamma
 Packing: 1pc/Sealed Aluminum Foil Bag, 100pcs/box, 5000pcs/ctn



FY0903 Surgical Handle

Size: 3#, 4#
 Index stainless steel metal handles are supplied non-sterile and individually wrapped
 Packing: 1pc/PE Bag, 10pcs/box, 500pcs/ctn



FY0904 Surgical Scalpel

Standard sheath guarded and blade protected sterile disposable scalpels offer convenience in a number of user environments: on the ward, in emergency or trauma units or even out in the field with paramedics
 Material: Carbon Steel or Stainless Steel
 Sterilization: Gamma
 Packing: 10pcs/Box, 500pcs/ctn



FY0905 Steel Lancet

Low cost: simple, flat, steel lancet
 Sample volume: excellent blood flow
 High quality: hardened surgical steel
 Minimum discomfort: modern asymmetric design
 Sterile: sterilized by gamma radiation
 Size: Adult 55mm
 Packing: 5pcs of packaged, 200pcs/Box, 20000pcs/ctn



FY0906 Twist Lancet

Conventional twist-off cap design
Sterilized by gamma-radiation
Smooth tri-bevel needle tip for comfortable sampling experience
Versatile with most lancing devices
Size: 21G, 23G, 26G, 28G, 30G
Packing: 100pcs/Box, 20000pcs/ctn



**FY0908 Lancing Device (Releasable)
FY0909 Lancing Device**

The user can install the lancet with the device easily.
And the device will work with the lancet smoothly without any blockage. The lancet won't fall on its own when to fire the device. Retracting automatically after firing, the launcher is safe and reliable.



FY0910 Safety Lancet

Ergonomic design offers a more secure hold when to activate the device
Tri-bevel needle tip ensures the sharpness and accuracy
The spring buffer the puncture making the procedure a more comfortable experience
The needle retracts immediately after use
Sterilization : Gamma
Size: 21G, 23G, 26G, 28G, 30G



FY0911 Pressure Activated Safety Lancet

Needle is fully shielded before and after use
Pre-loaded, pressure activated design
Ensures a good blood sample every time
Ultra sharp tri-bevel needle or bi-bevel blade, high speed penetration
Variety of needle specifications to suit most capillary blood requirements
Sterilization : Gamma
Size: 21G, 23G, 26G, 28G, 30G



FY091201 Button Activated Safety Lancet

Needle tip is sealed with cap through injection moulding
Pre-loaded, button operated design
Single-use only, sterilised by gamma irradiation, cannot be reused
High quality, ultra sharp tri-bevel needle, high speed penetration
Variety of needle specifications to suit most capillary blood requirements
Sterilization : Gamma
Size: 21G, 23G, 26G, 28G, 30G



FY0923 Contact-Activated Safety Lancet

Safe: Self-destruct device to ensure the needle is well protected and hided before and after use
Painless: Unique trigger will press the nerve ending, which can reduce the subject's feeling from the puncture

Cat. NO.	Needle Shape	Puncture Depth	Color
FY092301	21G(0.8mm)	1.8mm	Dark Green
FY092302	21G(0.8mm)	2.2mm	Orange
FY092303	21G(0.8mm)	2.8mm	Pink
FY092304	23G(0.6mm)	1.8mm	Light Green
FY092305	26G(0.45mm)	1.8mm	Yellow
FY092306	28G(0.36mm)	1.8mm	Purple
FY092307	30G(0.3mm)	1.8mm	Dark Gray

FY0913 / FY0914 / FY0915 Medical Razor

• Single or double blades • All razors come with the highest quality stainless steel blade or carbon steel, Comfort-coated for smooth shaves



FY0913



FY0914

FY0915



FY1501 Thermal Emergency Blanket

Material: PET Film
Waterproof, use to offset hypothermic reactions, reflective surface provides high visibility, reflects sun and heat for cooling in hot climates, space-age design retains up to 90% of body heat, Assists in prevention of post-accident shock and first aid application
Size: 132x210cm, 140x220cm, 160x210cm
Color: Gold, Sliver, Gold+Sliver



FY1502 Aluminum Mouldable Splint

Aluminum plate Macromolecule foam material Nontoxic and nonbenzenoid
Free of effect of temperature and climate Type available: Fold, Unfold and Roll
Normal size: 5cmx12cm, 5cmx15cm, 7cmx11cm, 11cmx44cm, 11cmx91.5cm, 11cmx100cm

FY2209 Disposable Electrosurgical Pencil



FY220901

Button Shape: Round
Hand control
Cable Length: 3m
Exposed electrode length: 40mm



FY220902

Button Shape: Oval
Handle Length: 170mm



FY220903

Foot control
Cable Length: 3m
Exposed electrode length: 40mm



SIZE:5x5cm

SIZE:2.5x5cm

FY2219 Electrosurgical Pencil Tip Cleaner




FY2220 Electrosurgical Pencil Holster

FY220904 Smoke Evacuation Electrosurgical Pencil

FY220905 Disposable Rocker Switch Pencil

Cable Length:3m

Teflon Coated (non-sticking coating prevent tissue from adhesive)



Disposable Electrosurgical Pencil (Tips)

Blade 16cm

Blade 14cm

Blade 12cm

Blade 10cm

Blade 9cm

Blade 7cm

Ball 8cm ϕ 6mm

Needle 7cm

Needle 15cm

Teflon Coated Available

Surgical Suture




Thread Length: 45cm, 60cm, 75cm, 90cm, 150cm etc.
Thread Diameter: 8/0,7/0,6/0,5/0,4/0,3/0,2/0,1/0,1,2,3 etc.
Needle Length: 6mm,8mm,12mm,18mm,35mm,40mm,50mm etc.

Needle Cross Section

Needle Shape	Point Type	Symbol
	Round Bodied	
	Curved Cutting	
	Reverse Cutting	
	Taper	
	Blunt Point	
	Precision Point-Reverse Cutting	
	Micropoint Curved Spatula	

Curvature



Suture Materials

Absorbable
FY1001 (CC)



Chromic Catgut

FY1002 (CP)




Plain Catgut

FY1007 (PGA)



Polyglycolic Acid


FY1008 (RPGA)



Polyglactin 910


Non-absorbable

FY1003 (SK)



Silk (Braided)

FY1004 (NL)



Nylon Monofilament

FY1005 (PB)



Polyester (Braided)

FY1006 (PM)



Polypropilene



EXAMINATION



FY0802 Surgical Gloves

Size: 6.0#, 6.5#, 7.0#, 7.5#, 8.0#, 8.5#, 9.0#
 Powder/Powder free
 1pair/pouch, gamma sterile
 AQL: 1.5
 Packing: 50pairs/box, 500 pairs/ctn



FY0803 Latex Examination Gloves

Size: S,M,L
 Powder/Powder free
 100pcs/box, non-sterile
 AQL: 2.5
 Packing: 100pcs/box, 1000 pcs/ctn



FY0804 Nitrile Examination Gloves

Size: S,M,L
 Powder free
 100pcs/box, non-sterile
 AQL: 1.5/2.5
 Packing: 100pcs/box, 1000 pcs/ctn



FY0805 Vinyl Gloves

Size: S (3.0g,3.5g,4.0g) ; M (3.5g, 4.0g, 4.5g) ;
 L (4.0g,4.5g,5.0g)
 Powder free
 100pcs/box, non-sterile
 AQL: 4.0
 Packing: 100pcs/box, 1000 pcs/ctn
 Color: Blue, White



FY0806 PE Gloves

Size: S,M,L
 Thickness: 0.018mm, 0.02mm, 0.022mm etc.
 Material: HDPE/LDPE
 Packing: 100pcs/box, 10000 pcs/ctn




FY0809 TPE Gloves

TPE gloves have stretch resistance and soft texture design
 Thickness: 0.20mm-0.45mm.
 Weight: 1.6-3.5 g/pc.
 Size: S, M, L, XL.
 Color: transparent white, blue etc.
 Applied to hospitals, restaurants, hotels, shops and ordinary household etc.

FY2201 Vaginal Speculum

The vaginal speculums are made from clear, non-toxic polystyrene plastic. The edges of the speculum are smooth and rounded, and very comfortable for patients. Also, they are easy to insert and manipulate with one hand only.
 Vaginal Speculum with Smoke Tube also available.
 Individual packaging: OPP, PE, Blister, Paper-foil;
 The speculums are available in different sizes, L, M, S, Even XS for some types.



FY220101 Side Screw Type



FY220102 Fastener Type



FY220103 Spanish Type



FY220104 Middle Screw Type



FY220105 French Type



FY220106 Push Pull Type



FY221202 Disposable Anoscope with Light


Disposable anoscope with integrated LED light source
 Single use, cable-less and self-illuminating proctoscope designed to optimize mobility during diagnosis and treatment.
 Button Battery SR-LR41, 1.5V
 Packing: Individually sterilized packing



FY220201 Cyto Brush



FY220203 Disposable Sampling Swab



FY220205 Disposable Cervical Brush



FY220206 Cervical Scraper



FY220207 Cervical Scraper



FY220208 Cervical Cell Sampling Spoon



**FY2101
Wooden Tongue Depressor**

Size:Adult (150mm × 18mm × 1.6mm)
Child (140mm × 14mm ×1.6mm)
EO Sterile / Non-Sterile
Packing: Single medical paper wrapped/Bulk
100pcs/box, 5000pcs/ctn



**FY2102
Cervical Scraper (Wooden)**

EO Sterile / Non-Sterile
Packing : Single medical paper wrapped/Bulk
100pcs/box, 5000pcs/ctn



**FY1901
HCG Test Strip**

Size: 2.5mm, 3.0mm, etc
Response time: 5 minutes
Sensitivity: 25mIU/ml
Specificity: 99%
Packing: 1pc/aluminum foil bag



**FY1902
HCG Test Cassette**

Size : 3.0mm, 3.5mm, etc
Response time: 5 minutes
Sensitivity : 25mIU/ml
Specificity : 99%
Packing : 1pc/aluminium foil bag



**FY1903
HCG Midstream**

Size : 4.0mm, 4.5mm, 5.0mm, etc
Response time: 5 minutes
Sensitivity : 25mIU/ml
Specificity : 99%
Packing : 1pc/aluminium foil bag



**FY2205
Gynecological Kit**

Configuration :
1pc Speculum
1pc Gynaecology Pad(40*60cm)
1pair Latex Gloves
1pc Cervical Spatula



FY1904 Condom

Specification:
Surface: smooth, dotted, ribbed lubricated
Size: L, M, S (55mm±2、 52mm±2、 49mm±2).
Be available with spermicide nonoxonyl 9
Be available with various fruit smell,such as apple/ banana/ strawberry/ cherry etc.
Be available with square or rectangle packing ;
If quantity is large enough, the condom is available with various color such as natural/ red/pink/green. etc.
Specification is available with EN600, ISO4074, ASTM, WHO, SABS, etc.



**FY1907
Ultrasound Probe Cover
(Telescope Folded&3D Dimensional)**

Product details:
Introduction:
It is made of PU material, apply to ultrasound probe and other surgical instruments. It could protect both medical instrument and patient.
Components:
Telescopic folded cover 1 pc
PE tapes 2 pcs
Rubber bands 2 pcs
Ultrasonic gel 1 pc
Feature:
1. PU material, dust-proof, water-proof, disposable and sterile
2. To prevent the equipments from contaminated and infected in surgery
3. Reduce the risk of cross-infection



MF207 PAP Smear Test Kit

The pap smear kit is used for efficient and adequate collection of uterine cervical specimens .
The enclosure package promotes specimen safety during collection and transportation.
Contents: Cervical brush, Cotton applicator, Wooden cervical scraper, Glass slide, Hard Blister Packing
1pc Glass Slide
1pc Cyto/Cervical Brush
1pc Wooden Applicator, 15cm
1pc Wooden Cervical Scraper
Hard Blister Packing

**FY9001
Lubricating Jelly Series**

The human body lubricant is composed of glycerin, polyethylene glycol, carbomer, methyl paraben, propyl 4-hydroxybenzoate, sodium hydroxide, and distilled water.
Suitable for lubrication of dry medical devices entering the human cavity.

Cat. NO.	Gram	Description	Qty/box	Qty/Ctn
FY900101	3g	0.8% glycerol	144/pcs	12/pcs
FY900102	57g	0.8% glycerol	12/pcs	6/pcs
FY900103	113g	0.8% glycerol	12/pcs	6/pcs
FY900104	5g	0.8% glycerol	144/pcs	4/pcs




WOUND DRESSING



FY1101 Gauze Swabs

The pads are folded all by machine, pure 100% cotton yarn ensure the product soft and adherent. Superior absorbency makes the pads perfect for absorbing blood and exudates.
In accordance with the customers' requirements we can produce different kinds of pads, such as folded and unfolded, with x-ray and without x-ray. The adherent pads are perfect for operation.
General specification: 2"×2", 3"×3", 4"×4"
PLY: 8ply, 12ply, 16ply
Mesh: 19×10, 19×15, 20×12, 26×18
EO Sterile / Non-Sterile



FY1105 Non-woven Swabs

Non-woven swabs from 70% polyethylene and 30% terylene.
The pad is designed to wick away fluids and disperse them evenly.
Product have been cut like "O" and "Y", so it is easy to use.
It is mainly used to absorbing blood and exudates during operation and cleaning the wounds.
General specification: 30g/m², 40g/m²;
2"×2", 3"×3", 4"×4"
PLY: 4ply, 8ply
EO Sterile / Non-Sterile



FY1103 Laparotomy Sponge

Laparotomy Sponge for absorbing blood and exudates or for cleaning wound.
Cotton Yarn: 40's x 40's
Type: Prewashed or non-prewashed.
It can prevent the wound from being infected. Our laparotomy sponges are sewn and cross-stitched.
General specification: 45×45cm/30×30cm,
PLY: 8ply, 16ply
Mesh: 19×15/26×18
EO Sterile / Non-Sterile



FY1102 Gauze Roll

Material : 100% cotton
Cotton Yarn: 40's×40's
Weaving according to customer requirements, high pressure, high temperature degreasing, bleaching, refined but become, with x-ray or without x-ray
General specification: 36"×100yard-4ply, 36"×100m-4ply, 36"×1000m-4ply;
Mesh: 19×11/19×15



FY1104 Gauze Bandage

Material : 100% cotton
 Cotton Yarn: 40's×40's
 Weaving according to customer requirements,
 high pressure, high temperature degreasing,
 bleaching, refined.
 Packing : 1roll/white paper, 1 doz
 blue kraft paper packing
 General specification:
 2" ,3" ,4" ,6" ×4.5m/5m/10m;
 Mesh:19×15,20×12,26×18
 EO Sterile / Non-Sterile



FY1106 Cotton Roll

Absorbent cotton is pure 100% cotton fabric,
 degreased and bleached by advanced way.
 The cotton is perfect for beauty; light drainae
 wounds, the alcohol sterile, cleaning the medical
 instruments.
 Packing become rolls by white and blue paper
 or PE bag.
 General specification:50g,100g,200g,250g,454g,
 500g,1000g
 EO Sterile / Non-Sterile



FY1109 Gauze Ball

100% pure cotton yarn of our product ensures
 a superior absorbency.
 Types with x-ray or without x-ray.
 It is mainly used to absorbing blood and exudates
 during operation and cleaning the wounds.
 General specification: Dia.:3cm, 4cm, 5cm;
 Mesh:19×15, 26×18
 EO Sterile / Non-Sterile



FY1107 Zig-zag Cotton

Absorbent cotton is pure 100% cotton fabric,
 degreased and bleached by advanced way.
 Packing: sealed poly bag, drawstring bag
 or self-sealing bag
 General specification: 25g,50g,60g,100g,250g,500g
 EO Sterile / Non-Sterile

FY1108 Cotton Ball

Soft & Pure balls are also highly absorbent, offering
 outstanding cleansing properties as they absorb
 fluids extremely quickly.
 Made from 100% pure cotton
 Gentle to the skin
 General specification: 0.2g, 0.3g, 0.5g, 1.0g
 EO Sterile / Non-Sterile



FY1110 Cotton Applicator(Cotton Swab)

Size: 3", 4", 6", 8"
 Material: Wooden or Plastic Sticks with 100% absorbent cotton
 EO Sterile/Non-Sterile



FY1318 Abdominal Pads

The pads consist of a highly absorbent inner (cotton), enclosed
 in a soft and comfortable non-woven fabric or Gauze.
 They are designed to have high absorbency, which makes
 them ideal for managing large amounts of wound exudate
 or fluid drainage. They can effectively wick away moisture,
 keeping the wound area clean and dry.
 With/Without detectable line
 Size: 5"x9", 5cmx5cm,10cmx10cm, 10cmx20cm,
 20cmx20cm,20cmx30cm,etc



FY1412 Combined Dressing Roll (Gamgee Roll)

Roll dressings consist of a highly absorbent inner(cotton),
 enclosed in a soft and comfortable non-woven fabric or Gauze.
 The roll can be cut to sizes that suit.
 With detectable line or without detectable line.
 Different sizes are available.



With Detectable Line

Without Detectable Line

FY1401 Dental Cotton Roll

The dental cotton roll is made of 100% cotton, It is widely used on dentistry and dental cleaning.
 Size: 0.4g, 0.5g, 0.7g
 Size: 8 × 38mm, 10 × 38mm, 12 × 38mm
 EO Sterile / Non-Sterile



FY1406 Oral Cleaning Swab

Material: PP and sponge
 Color: green, red, blue, white or others;
 Product Specification: 10cm, 15cm, 18cm, or as your request.
 Gently cleans, moistens and refreshes mouth, teeth, gums and disinfects the wounds.
 Raised foam ridges remove particles and mucous for enhanced oral hygiene



FY1402 Saliva Ejector

Material: PVC
 Length: 15cm
 Packing: 100pcs/bag, 20bags/ctn,



FY1419 Disposable Dental Air Water Syringe Tips

Directly interchangeable with most standard style metal air/water Syringe tips
 No sharp edges for patient comfort
 Unique quad air core design to ensure 100% dry air
 Different color is available



FY1404 Dental Bib

Material: 1ply or 2ply 17g/Sqm Tissue + 1ply 13g/Sqm Poly Film
 Color: White/green/pink/blue, etc.
 Packing: 125pcs/bag, 4bags/ctn



FY1414 Bib Clip

Cat.No.	Name	Color	Material
FY141401	Bib Clip (all metal)	/	Made of stainless steel bead chain and clip
FY141402	Bib Clip (plastic spring)	Red, Blue	Made of plastic spring and stainless steel clip Good elastic suit for different patients
FY141403	Autoclavable Dental Bib Holder	Blue, Green, Light Blue, Purple	Made of medical grade silicone body and stainless steel clip Fully autoclavable up to 135°C



FY1408 Disposable Plastic Forceps

Material: Plastic
 Color: red, blue, green, yellow, white or as customers' request.
 Dimension: 10.8cm, 11.5cm, 12.5cm, 13 cm, 14cm, 18 cm, etc.



FY1409 Plastic Haemostatic Forceps

Material: Plastic/ABS
 Color: red, blue, green, yellow, white or as customers' request.
 Dimension: 14×7.5cm



FY1201 Spandex Elastic Bandage

Material: cotton and spandex
 Color: natural, flesh or bleach
 Weight: 60g/m², 70g/m², 75g/m², etc.
 Width: 5cm, 7.5cm, 10cm, 15cm
 Length: 4m, 4.5m, 5m; etc



FY1202 PBT Elastic Bandage

Material: cotton and polyester
 Color: natural or bleach
 Weight: 25g/m², 28g/m², 30g/m²
 Width: 5cm, 7.5cm, 10cm, 15cm
 Length: 4m, 4.5m, 5m; etc



**FY1203
Cotton Elastic Bandage**

Material: 100% cotton
Color: natural or bleach
Weight: 60g/m², 70g/m², 75g/m², etc.
Width: 5cm, 7.5cm, 10cm, 15cm
Length: 4m, 4.5m, 5m; etc



**FY1204
High Elastic Bandage**

Material: rubber and polyester
Color: white, flesh or bleach
Weight: 80g/m², 100g/m², etc.
Width: 5cm, 7.5cm, 10cm, 15cm
Length: 4m, 4.5m, 5m; etc



**FY1205
Elastic Net Bandage**

Material: rubber and polyimide
Color: white
Weight: 120g/m², etc.
Size: 0.8cm×25m, 1.7cm×25m, 2cm×25m,
2.2cm×25m, 2.6cm×25m,
3.8cm×25m, 4.4cm×25m, 5.5cm×25m,
6.8cm×25m, 8cm×25m



**FY1206
Tubular Bandage**

Material: rubber and cotton
Color: natural
Width: 4.5cm, 6.25cm, 6.5cm, 6.75cm, 7.5cm,
8.75cm, 10cm, 12cm, 17.5cm, 21.5cm, 32.5cm
Length: 10m



**FY1220
Stockinette Bandage**

Material: 100% Cotton
Color: Natural white, or other color color customer's demand
Width: 5cm, 6cm, 7.5cm, 8cm, 10cm, 15cm, 20cm
Length: 20m

**FY1207
Non-woven Self-adhesive Bandage**

Material: non-woven
Weight: 35g/m², 40g/m², 45g/m², etc.
Color: 12 colors
Width: 5cm, 7.5cm, 10cm, 15cm
Length: 4m, 4.5m, 5m; etc



**FY1208
Cotton Self-adhesive Bandage**

Material: cotton
Weight: 75g/m², 80g/m², 85g/m², etc.
Color: 12 colors
Width: 5cm, 7.5cm, 10cm, 15cm
Length: 4m, 4.5m, 5m; etc



**FY1209
Plaster of Paris Bandage**

Material: cotton and plaster
Weight: 390g/m², 410g/m², etc.
Color: white
Width: 5cm, 7.5cm, 10cm, 15cm
Length: 2.7m, 3.6m, 4.5m, etc



**FY1217
Orthopedic Casting Bandage**

Description: Fiber glass casting tape consists of a lightweight knitted polyester fabric, impregnated with a water receptive polyurethane resin.
Setting time: 3-5 minutes
Weight-bearing: 20 minutes
Excellent conformability & mouldability
Tack free 100% transparent to X-rays
Very good water vapour permeability
With elasticity, and flexibility
Blue, Green, White, Pink, Red color
Width: 5cm, 7.5cm, 10cm
Length: 2.7m, 3.6m, 4.5m, etc



**FY1218
Cast Padding**

Material: natural cotton, polyester, etc
Color: white
Width: 5cm, 7.5cm, 10cm, 15cm
Length: 4m, 4.5m, 5m, etc





**FY1211
Zinc Oxide Plaster**

Made by using cotton cloth as a base material and coating adhesive cement mainly made of zinc oxide and rubber or medical hot-melt adhesive.
Packing: simple package, plastic can, galvanized can
The product is high in stickiness, and can be cut into a saw-tooth shape
Color: white, flesh
Width: 1.25cm, 2.5cm, 5cm, 7.5cm, 10cm, etc.
Length: 4.5m, 5m, 10yard, 10m, etc.



**FY1212
Non-woven Surgical Tape**

The non-woven surgical tape mainly uses low-sensitization medical hot-melt adhesive, the back face is coated with an uncoiling agent, and the back lining material is mainly made of soft, light, thin and air-permeated nonwoven cloth serving as a base material.
Packing: simple package, plastic can, with cutter
Color: white, flesh
Width: 1.25cm, 2.5cm, 5cm, 7.5cm, 10cm, etc.
Length: 4.5m, 5m, 10yard, 10m, etc.




**FY1213
PE Surgical Tape**

The PE surgical tape mainly uses low-sensitization medical hot-melt adhesive, the back face is coated with an uncoiling agent, and the back lining material is mainly made of soft PE films serving as a base material.
Packing: simple package, plastic can, with cutter
Color: transparent
Width: 1.25cm, 2.5cm, 5cm, 7.5cm, 10cm, etc.
Length: 4.5m, 5m, 10yard, 10m, etc.



**FY1214
Silk Surgical Tape**

The silk surgical tape mainly uses low-sensitization medical hot-melt adhesive, the back face is coated with an uncoiling agent, the back lining material is mainly made of soft, light and thin silk cloth serving as a base material, and the adhesive is coated in a shape of mesh.
Packing: simple package, plastic can, with cutter
Color: white
Width: 1.25cm, 2.5cm, 5cm, 7.5cm, 10cm, etc.
Length: 4.5m, 5m, 10yard, 10m, etc.



**FY1222
Kind Removal Silicone Tape**

Security: Offers reliable and pliable fixation, remains in place until you decide otherwise
Comfort: Removes cleanly, without disrupting fragile skin layers or causing patients undue pain
Ease of Use: Can be repositioned and neatly torn by hand
Economical: Offered in single patient roll sizes
Width: 1", 2"
Length: 1.5yard, 5.5yard

**FY1215
Sports Tape**

The sports tape uses cotton cloth serving as a base material and adhesive cement mainly made of zinc oxide and rubber, and the adhesive cement is coated in a shape of strip.
Different colors available
Width: 2.5cm, 3.8cm, 5cm, etc.
Length: 10yard, 15yard, etc.



FY1225 Kinesiology Tape

Material: Waterproof Elastic Fabric
Longitudinal extension rate: 1:1.5~1:1.8
Coating: Curve Adhesive Coating
Size: Width: 2"
Length: 16.4 feet, 13feet
Color: Black, Light Blue, Green, Orange, Pink, Purple, Red, Dark Blue, Skin, Yellow, etc



**FY1226
Elastic Adhesive Bandage**

Heavy-adhesive elastic bandage has good stickiness, breathable, comfort, no degumming and residue.
Its cotton material is soft and highly applicable and its tensile strength can be freely adjusted, yellow middle line design.
Used for sprains and strains; fixing the cold or hot compress bag; controlling edema after varicose veins; compression hemostasis; treatment of animal's disease.
With film or without film is available.
Width: 5cm, 7.5cm, 10cm
Length: 4.5m



**FY1216
Cold Wrap Elastic Bandage**

Cold Wrap Elastic Bandage is a cold bandage which uses a revolutionary new technique for treating injuries, muscle pain and strains. This bandage reduces recovery time from injuries by combining both the cooling and compressing process in one easy to use step.
Color: brown, blue
Size: 7.5cm x 3.2m, 10cm x 3.2m



**FY1301
Drilled Plaster**

The plaster is perforated zinc oxide plaster, using special processing to increase the viscosity and permeability of products, especially suitable for all kinds of environment.
It can be cut to required size by scissor, comfortable and breathable, zinc oxide glue provides strong fixation.
Width: 1.25cm, 2.5cm, 5.0cm, 7.5cm, 10cm.
Length: 5m, 10m





**FY1302
Healing Capsicum Plaster**

The product is mainly formed by two parts: a base material (cotton material) and separate paper. The product has the features of air permeation, moisture permeation, tight fixing, High adaptability, convenient use, etc



**FY1316
Catheter Securement Device**

The Catheter Securement Devices are aimed to improve and protect patient safety. They can secure the I.V. Catheter, Foley Catheters, Nasal Gastric tubes and Various d of Catheters & Tubes tightly, help to reduce the risk of catheter dislodgement. Soft and flexible fabric design improves patient comfort. Provides superior securement for both horizontal and vertical lifting accidental pulls. Simple to apply, inspect and adjust. Hypoallergenic, breathable and latex-free to reduce the risk of allergic reactions and skin irritation.



**FY1303
Wound Plaster**

Provides wound care for minor cuts, scrapes & burns
Size: 72x19mm, 70x18mm, 76x25mm, Dia22mm, Dia25mm etc
Material: PE, PVC, EVA, PU, Fabric, Cotton, Etc



**FY1304
Eye Pad**

There are two kind of Eye pad. One is non adhesive 100% Cotton eye pad. The other is made of an elastic non-woven dressing with absorbent pad. Both dressing are Specially indicated for strabismic, weak sighted therapy and nursing after ocular operation. Color: White, Skin

**FY1633
Non-woven Dressing Roll**

Material: Spun-laced Fabric
Composed: Non-woven Fabric + Medical Pressure-Sensitive Glue + Release Paper.
Release Paper Cutting Type: Straight Cutting/Curved
Package: Individual Box
Size: Width: 1", 2", 3", 4", 6", 8", 10"
Length: 5m, 10m, 20m
Size can be customized according to customer's requirements.



**FY1642
PU Dressing Roll**

Material: PU (Transparent, Conformable, Unique roll format OPSITE Film)
Composed: PET packing + PU Film + Medical Pressure-Sensitive Glue + Release Paper.
The film is waterproof and aids in the prevention of contamination.
Release Paper Cutting Type: Straight Cutting/Curved
Package: Individual Box
Size: Width: 1", 2", 3", 4", 6", 8", 10"
Length: 5m, 10m, 20m
Size can be customized according to customer's requirements.



**FY1305
Non-woven Adhesive Pad**

The external layer of the dressing pad is formed by spun-laced nonwoven cloth, the internal layer is formed by needle-punched nonwoven cloth, and the surface is formed by release paper; the base material is soft, and the product should be sterilized and is mainly used for wound dressing in the department of surgery of a hospital.
Size: 6cm×8cm, 10cm×10cm, 10cm×20cm, etc.



**FY1317
PU Adhesive Dressing**

PE square backing design
Capable of preventing tension on the skin during the operation, and reducing the probability of skin injury.
Super high breathability, waterproof and skin-adaptive
The color of the imported ultrathin frosty polyurethane (PU) with super high breathability is close to that of the skin, transparent, breathable and waterproof. It can be fully attached on the skin and convenient for observation of the injection position and wound.
Size: 6cm×8cm, 10cm×10cm, 10cm×20cm, etc.





**FY3201
PU IV Dressing Pad**

The polyurethane (PU) film is completely transparent and Permeable, making it easy to observe the injection site and wound conditions. The unique U-shaped incision design facilitates the fixation of indwelling needle in butterfly wing shape or with medicine pot. Ultra-thin design for perfect fitting on the skin. Special ear shape makes paper frame stripping easier. Additional record and breathable tape for fixation make it easy for operator to record patient name, time and bed No., and can be used to strengthen the fixation for quality management.
Size: 6cm x 7cm, 10cm x 12cm, etc.



**FY1221
Wound Skin Closure**

Suitable for securing, closing and supporting small cuts and wounds
Material: Special non-woven fabric (with a line or not with a line)
Size: 3x75mm/6x75mm/6x100mm/12x100mm



**FY1325
Hydrocolloid Dressing With Border**

It is a kind of elastic dressing made of pressure sensitive adhesive with CMC(Carboxymethylcellulose) hydrophilic granules. They provide a moist and insulating healing environment which protect uninfected wounds while allowing the body's own enzymes to help heal wounds. accelerate wound healing process, effectively reduce pains and wound care frequency. Easy to move off, and no adhesive with wound, so effectively reduce the possible damage and pains.
Size: 15cm x 20cm, 20cmx20cm, etc



**FY1327
Trachea Foam Dressing**

Sterile, hydrophilic foam dressing with a circular aperture designed to fit neatly around a tracheostomy tube or other drain or stoma; polyurethane foam wound contact surface has high absorption capacity and a vapor permeable, water- and bacteria-resistant polyurethane film outer layer.
Size: 9cmx9cm, 10cmx10cm, etc

**FY1328
Silicone Foam Dressing with Border**

It is used in surgical wounds with more exudate and urology, plastic surgery, orthopedics, and ophthalmology. Silicone gel foam dressings have a five-layer structure:
1. Silicone gel layer
2. Foam absorbent layer
3. One-way flow layer
4. Highly absorbent liquid storage layer
5. Semi-permeable membrane layer
Size: 7.5cm x 7.5cm, 10cm x 10cm, 15cm x 15cm, 10cm x 20cm, etc.



**FY1329
Silicone Foam Dressing**

It is suitable for wounds with moderate to a large amount of exudate, reducing trauma, reducing pain, absorbing a sufficient amount, and promoting wound healing.
Features:
1. Minimize wound pain.
2. Minimize damage to the surrounding skin.
3. Prevent wound maceration.
4. Keep the wound moist to speed up wound healing.
Size: 5cm x 5cm, 7.5cm x 7.5cm, 10cm x 10cm, 10cm x 20cm, etc.



**FY1310
Vaseline Gauze Swab(Paraffin Gauze)**

The product is made from medical degreased gauze or non-woven together with Vaseline. It can lubricate the skin and protect the skin from cracks. It is widely used on clinic.
Packing: Individual packing or 36 pcs per plastic can



**FY130601 / FY130602
Alcohol Prep Pad(Alcohol Swab)**

Contents: 70% Isopropyl Alcohol
Non-Woven Size : 65x30mm, 65x58mm, etc
Type: 2ply, 4ply, etc.
Packet Size: 50x50mm



**FY1314
Adhesive Remover Pad**

Helps to remove adhesive tape, adhesive residue and pouches from skin
Size: 65x30mm, etc
Packing: individual Packet, 50pcs per box





FY1307
Povidone-Iodine Prep Pad

Contents: 10 % Povidone-Iodine Solution
Non-Woven Size : 65×30mm,65×58mm,etc
Type: 2ply,4ply,etc.
Packet Size: 50×50mm



FY1308
Sting Relief Medicated Pad

Contents: Antipruritic solution
Non-Woven Size : 65×30mm,65×58mm,etc
Type: 2ply,4ply,etc.
Packet Size: 50×50mm



FY1309
Povidone-Iodine Swabstick

Contents: 10 % Povidone-Iodine Solution
Swabstick Size: 70mm,100mm
Packet Size : 110×40mm,140×40 mm
Package: 1pc or 2pcs or 3pcs per Airtight Sacket



FY1311
Alcohol Swabstick with Chlohexidine

Alcohol Swabstick Packets are useful for topical cleansing and skin prep prior to injections or venipuncture.
Contents: 70% Isopropyl Alcohol or 70% Isopropyl Alcohol + 2% Chlorhexidine
Swabstick Size: 70mm, 100mm
Packet Size: 110x40mm, 140x40mm
Package: 1pc or 2pcs or 3pcs per Airtight Sacket



FY1312 Wet Wipes

Specifications:
Material: Spunlace Nonwoven, ES hot air through Nonwoven, Thermal-bonded Nonwoven, Air-laid paper, Wet Strength Paper, Bamboo Fiber.
Weight: 25-80gsm for your choice
Size: 15x20cm, 16x20cm, 18x20cm, 20 x 30cm, etc.
Wipes count: 1 to 200pcs for your choice.
Folding Size: Z fold, M fold.
Formula: Gentle and safe, Alcohol Free, all kinds of fragrance for your choice.
Fragrance: Lemon, Lavender, Camellia, etc.
Function extracts for your choice, such as Aloe Vera & Vitamin E and so on.
Packaging type: PE/PET film, plastic canister/tube, or plastic box, etc.

**MEDICAL
NON-WOVEN
& ACCESSORIES**



**FY1601
Round Nurse Cap**

Material: Non-woven PP
Weight: 10g/m²,12g/m²,14g/m²,16g/m²,etc
Size: 18", 21", etc
Color: white, blue, green, etc



**FY1602
Clip Cap**

Material: Non-woven PP
Weight: 10g/m²,12g/m²,14g/m²,etc
Size: 20", 21", 22", etc
Type: single elastic/double elastic
Color: white, blue, green, etc



**FY1603
Doctor Cap**

Material: Non-woven PP
Weight: 20g/m², 25g/m², 30g/m², etc
Size: 61×14cm,64×23cm, etc
Type: tie/elastic back
Color: white, blue, green, etc



**FY1604
Net Cap**

Made from high quality nylon net fabric material
Weight: 2.7g/m², etc
Size: 19", etc
Color: white, blue, green, etc



**FY1609
Surgical Hood**

Material: Non-woven PP
Weight: 16g/m², 20g/m²,25g/m²,30g/m²,etc
Size: 30×124cm, etc
Type: tie
Color: white, blue, green, etc




**FY1605
Non-woven Face Mask**

Material: Non-woven PP or non-woven PP and filter paper
Specification: ear-loop/tie, 2ply/3ply/4ply
Size: 17.5×9.5cm
Color: white, blue, green, etc



**FY1632
Non-woven Face Mask with Charcoal**

Description:
Conforms to American NEL SON N95 test standards. Adjustable nose seal, easy to fold and carry. Mini-fiber high electrostatic filter, filtration efficiency of dust & powder is over 95%. Repel continually contamination of dangerous blood and saliva continually.
Specification: Ear-loop type: 4 ply
Detailed weight: 4 ply: 25g/m²+20g/m²+40g/m²+25g/m²
Size: 18×9cm
Color: white, green, blue etc.
Packing: 50pcs/box, 40boxes/carton



**FY1606
Paper Face Mask**

Material: high quality wet-strength paper
Specification: ear-loop and 1ply/2ply
Size: 20 × 7cm
Color: white



**FY1608
Face Mask with Shield**

Composition : non-woven face mask and PET shield (prevent mist)
Specification : ear-loop and 3ply/4ply
Size : 17.5x9.5cm
Color : white,blue,green,etc



**FY1629
Face Shield**

Material : PET
Type : elastic double-sided anti-fog
Pack : 1pc/pe bag



**FY1301
Goggles**

Impact resistant PC lens, Lens thickness: 2mm
Integrated side legs
High strength structure
Integrated wide side shields



**FY3501
Phototherapy Mask Eyecare**

Feature:
Easy-to-use design
Bigger eye pad for less pressure and better UV protection
Safe and comfortable
Special material to make the eyemask hard to be fluffy
Single patient could be used 3-5days during phototherapy treatment

The Eyecare is available in four sizes and come in packs of 20
N-A head dicumference 30-38cm
N-B head dicumference 24-33cm
N-C head dicumference 20-28cm
T head dicumference 24-35cm



**FY160702
N95 Mask (with/without valve)**

N95 is one of nine classifications for NIOSH certified particulate respirators.
"N" is for non-oil particulates;
"95" means the filtration efficiency will be no less than 95% when under the detecting conditions assigned by NIOSH.
Only those which meet the N95 standard and approved by NIOSH can be called NIOSH N95 masks.

Material: polypropylene (conforms to American NIOSH N95/FFP1/FFP2/FFP3 test standards)
Specification: ear-loop
Color: white



**FY160703
N95 Facemask (with/without valve)**

Material: polypropylene (conforms to American NIOSH N95/FFP1/FFP2/FFP3 test standards)
Specification: ear-loop
Size: without valve or with valve
Color: white



**FY160704
Folded Facemask (with/without valve)**

Material: polypropylene (conforms to American NIOSH N95/FFP1/FFP2/FFP3 test standards)
Specification: ear-loop
Size: without valve or with valve
Color: white

**FY1610
Surgical Gown**

Standard and Reinforced
Material: Spunlace, SMS, SU/US nonwoven and so on
Gram weight of material: from 35g/m² to 50g/m²
Color can be decided by customer
Whether sterilize is decided by customer
Packaging way can be bulk or individual
Size: S, M, L, XL, XXL, XXXL



**FY1611
Disposable Isolation Gown**

Material: SPP or PP+PE non-woven fabrics or as customers' requirement
Color: Blue, yellow, green and so on.
Size: S/M/L/XL/XXL
Details: knitted or elastic cuff or thumb type with velcro or ties on the neck and ties on the waist.



**FY1612
CPE Coat**

Material: CPE
Color: Blue
Size: 62g, 125×200cm, etc
Characteristics: prevent and isolate dust, particle, alcohol, blood Bacteria and virus invading



**FY1613
Lab Coat**

Material: non-woven PP/PE
Color: white/blue/yellow/green
Weight: 30g, 35g, 40g, 45g, etc
Size: M/L/XL/XXL



**FY1610R
Reinforced Surgeon Gown**

Specification:
Material: Spunbond Meltdown Spunbond (SMS)
Color: blue
Weight: 40gsm/45gsm/48gsm/50gsm/55gsm etc.
Size: S/M/L/XL/XXL
Spunbond meltdown spunbond(SMS) material gown breathable, quiet and supple,
High repellency, AAMI standard, sterile with comfortable design for operating theater
With knitted cuff, Velcro neck back closure, side waist band
Reinforced arm and additional liner at front
1 hand towel is available
1 wrap towel is available
Sterile method: non-sterile or EO sterile
Packing: 10pcs/poly bag or 1 pc/sterile bag





FY161501 Disposable Protective Coverall without Tape

Chemical protective coverall includes: liquid-tight protective clothing , spray-tight protective clothing protective clothing for use against solid particulates protective clothing against liquid chemicals and so on.
Material: PP/PP+PE/SMS and so on
Weight: 35-60g
Size: M/L/XL/XXL
Color: White
Application: provide reliable protective products for industry/lab/ nuclear radiation/ agriculture/ military combat readiness and so on



FY161502 Disposable Protective Coverall with Tape

Chemical protective coverall includes: liquid-tight protective clothing , spray-tight protective clothing protective clothing for use against solid particulates protective clothing against liquid chemicals and so on.
Material: PP non-woven fabric (30gsm) + Breathable film (30gsm)+glue(3gsm)
Size: M/L/XL/XXL
Color: White (blue adhesive tape)
Application: provide reliable protective products for industry/lab/ nuclear radiation/ agriculture/ military combat readiness and so on



FY1614 Disposable Scrub Suit

Size: S/M/L/XL/XXL and so on
Material: SS/PP/SMS 20g - 60g
Color: Blue/white/yellow/ dark blue
Packing: 1pc/PE bag, 50bags/carton



FY1622 PE Apron

Made of polyethylene, with adjustable tie
Color: white, blue, etc
Size: 110×70cm, 135×71cm, 160×90cm
Thickness: 0.018mm, 0.02mm and so on

FY1616 PE Sleeve

Made of PE, with comfortable elastic
Color: white, blue and red, etc
Thickness: 0.018mm
Size: 20×40cm, 22×42cm, 22×46cm



FY1617 PE Shoe Cover

Material: PE, with comfortable elastic opening
Color: blue, green, etc
Weight: 1.6g, 1.8g, 2.0g, 2.1g, 2.5g
Size: 36×15cm, 40×16cm
Thickness: 0.018mm, 0.02mm, 0.021mm and so on



FY1618 CPE Shoe Cover

CPE shoe cover has anti-slip particles on the surface, which is more anti-slip than PE shoe cover also stronger because it is thicker, and has the functions of heat-resistant oxygen aging, ozone aging, acid and alkali resistance, and corrosion resistance.
Color: blue, green, etc
Weight: 1.6g, 1.8g, 2.0g, 2.1g, 2.5g
Size: 36×15cm, 40×16cm
Thickness: 0.018mm, 0.02mm, 0.021mm and so on




FY1619 Non-woven Shoe Cover


Material: Non-woven PP
Color: white, blue
Regular/Non-skid
Weight: 30g/m², 40g/m², 45g/m², etc
Size: 36×15cm, 38×18cm, 40×16cm



FY1623 Non-woven Bed Sheet

Soft touch and harmless to your skin
100%biodegradable and environmentally friendly
SPEC: light blue , 60×180cm
30gsm other size and specification are also acceptable





**FY1624
Non-woven Bed Cover**

Disposable bed cover made of pp non-woven fabric helps you to be clean & hygiene on your travel and regular purpose
SPEC: 30g, 210×90×20cm, other size and specification are also acceptable



**FY1625
Non-woven Pillow Cover**

Disposable pillow covers are used widely in hotels hospitals, home, and airline and in all travels

Specification:
80gsm, 50×60cm, with 5cm gusset,
other size and specification are also acceptable
Packing: 100pcs/polybag, 1bag/ctn
Ctn size: 54×32×38 /100pcs



**FY1620
Under Pad**

Material: Non-woven fabrics+2ply medical paper + fluff pulp+PE film, non-woven fabrics+4ply medical paper+PE film, etc
with or without SAP
Size: 40×60cm, 60×90cm, etc
Weight can be decided by customer



**FY1621
Diaper**

First layer: Non-woven fiber, natural soft and breathable.
Second layer: full pulp+super absorbent polymer +classic tissue.
Breathable PE film with cloth like back Non-woven or casting film.
Standard Size:M: 650×800mm
L: 750×900mm, XL:840×1030mm

FY140702 Surgical Brush

Medical Brush has an effective, soft and skin-friendly brush head for the nails, knuckles and skin. Attached to the brush is a sponge designed for cleaning the hands and forearms.
In addition a separate nail cleaner is supplied in the packaging.
The Medica Brush saves on handling costs since brush, soap and nail cleaner are all supplied in one packaging.
In addition to a dry sterile version, one has a choice of several impregnated versions:
Medical Brush Dry Sterile
Medical Brush Chlorhexidine
Medical Brush Povidone Iodine




FY140701 Surgical Scrub Brush

Surgical scrub brush is used to significantly reduce the number of micro-organisms on the hands and forearms before performing surgeries. The nail cleaner feature makes removal of micro-organisms under the nail quick and easy
Dry Surgical Brush includes a variety of antimicrobial formulations, including CHG, PCMX, and Povidone Iodine
All kinds of spray cans are available



FY1630 Surgical Drape

Drapes can be with or without hole, it is available in different dimensions with different aperture sizes and shapes



FY2510 Basic Dressing Set

Dressing Tray	1Pc
Gauze swab 7.5cmx7.5cmx8ply	5Pcs
Tweezers	5Pcs
Hand Towel	3Pcs
Surgical Drap 50x60cm	1Pc
Limpet Bag	1Pc
Sterile Field 42x50cm	1Pc



FY2512 Dressing Set

Green plastic tweezers	1Pc
Non-woven swabs 30gsm, 7.5*7.5cm-4ply	10 Pcs
Sheet (paper+PE film) 38*45cm	1Pc
Hemostat Forceps	1Pc





- FY150301
Sterile Delivery Pack**
- Non-woven surgical gown 1pc
 - 86×80cm packing clothes 2ply (Non-woven and thin film) 1pc
 - 97×73cm sheets 3ply (Thin film, Non-woven and paper pulp) 1pc
 - 97×73cm Triangle trouser 3ply (Thin film, Non-woven and paper pulp) 2pcs
 - 66×43cm Non-woven cloth 2ply (Thin film and Non-woven) 2pcs
 - 66×43cm Non-woven cloth 1ply-1pc; Umbilical cord clamp 1pc
 - 8cm cotton sticker 2pcs
 - PE glove 1pair
 - Latex examination glove 1pair
 - 8×10cm Gauze swab 5pcs
 - Face mask 1pc
 - Cap 1pc



- FY150302
Sterile Delivery Pack**
- Instrument table cover 150×200CM 1Pc
 - Side drape 100×130CM 1Pc
 - Simple drape 100×130CM 1Pc
 - Hand towel 40×40CM 1Pc
 - Leggings drape 75×120CM 2Pcs
 - Clip 1Pc
 - Reinforced gown L 45gsm SMS 1Pc
 - Gauze 7.5×7.5cm 10 Pcs
 - Wrapping cloth 100×100CM 1Pc
 - Side drape 75×90CM 1Pc



- FY4201
Spinal Epidural CVC Sterile Pack**
- Fenestrated pe Drape 610 × 700mm, 1Pc
 - Luer Lock Syringe, 3ml and 5ml 2Pcs
 - Applicator Sponge, 3Pcs
 - Instrument table cover, 120×100cm 1Pc
 - Hand Towel, 30×40cm 2Pcs
 - Green Gauze 100×100mm 8Pcs
 - Filter Needle 19g, 2Pcs
 - Medical Plastic Tray, 24×25×4.5cm 1Pc
 - 1pack/blister packing, 30packs/ctn



- FY4202
Disposable Urethral Catheter Set**
- 10ml Lubricant Jelly 1Pc
 - Disposable Bag 1Pc
 - Dress Towel 1Pc
 - Disposable Forceps 1Pc
 - Disposable Container 1Pc
 - Fenestrated Drape 1Pc
 - Plain Drape 1Pc
 - Disposable Examination Gloves 1Pair
 - Large BPC Wool Balls 5Pcs
 - Tray 1Pc
 - Wrap 1Pc





**FY1701
Hot Water Bag**

Material: calcium carbonate≥30% rubber
Size: 500cc, 1000cc, 1750cc, 2000cc, etc
Packing: 1pc/PE bag




**FY1706
Ear Syringe**

Material: calcium carbonate≥25% rubber or PVC
Size: 30ml, 60ml, 90ml, etc
Packing: 1pc/PE bag



**FY170403
Hot/Cold Pack**

Material: PVC+carbamide + distilled water
Weight: 40g
Size: 9.8×8cm



**FY170401
Reusable Oval Ice Bag**

Material: Rubber
Max Length: 28cm
Max Width: 17cm



**FY170402
Hot/Cold Pack**

Material: PVC+carbamide + distilled water
Weight: 65g
Size: 10×10cm



**FY1705
Air Cushion**

Material: calcium carbonate≥30% rubber
Size: φ35cm, φ40cm, φ42.5cm, φ45cm, etc
Packing: 1pc/PE bag



**FY1710
Enema Rectal Syringe**

Material: Medical grade PVC&PE
Color: Any color available, customer customization

Cat. NO.	Volume	Length(mm)	Diameter	Weight
FY171001	41ml	89mm	47mm	24g
FY171002	89ml	137mm	60.5mm	44g
FY171003	143ml	153mm	70mm	66g
FY171004	224ml	168mm	82mm	102g
FY171005	275ml	193mm	86mm	142g
FY171006	347ml	205mm	90mm	167g

**FY1702
Fabric Ice Bag**

Made of canvas
Size: 6", 9", 11", etc
Packing: 1pc/PE bag



**FY1703
Instant Ice Bag**

Material: PE+PA infilling (carbamide + distilled water)
Weight: 160g, 280g
Size: 19.5×12cm, 24×12cm



**FY170404
Reusable Hot/Cold Pack**

Material: nylon taffeta infilling (glycerol)
Weight: 200g
Size: 20.5×15.5cm

**FY0233S
Silicone Enema Bag**

Material: Edible grade silicone
Type: Shiny surface or Mat surface
Size: 2L, 3L, 4L



**FY0702
Pressure Infusion Bag**

Use for emergency transfusion and invasive arterial pressure monitoring
Size: 500cc, 1000cc, 3000cc
Cuff: Disposable, Reusable
Sleeve Style: Nylon Mesh or Transparent TPU
Pressure Gauge: Piston gauze, Aneroid gauze, Waterproof piston gauze
Pressure Control Valve : Stopcock, Thumbwheel valve



FY2001
Feeding Bottle

Material: PC/Glass+Liquid silicone nipple
Description: 150ml, 225ml, 250ml, 280ml



FY2002
Umbilical Cord Clamp

Material: ABS, PE, PP
Color: blue, white, etc
Size: 5cm, 5.5cm



FY2301 Heat-Sealing Gusseted (Reel)

Size: 75mm(3")×100m, 100mm(4")×100m
150mm(6")×100m, 200mm(8")×100m
250mm(10")×100m, 300mm(12")×100m
350mm(14")×100m, 400mm(16")×100m



FY2302 Self-sealing Flat Pouch

Size: 50×130mm, 90×230mm
135×255mm, 190×330mm
250×370mm



FY2003
Umbilical Cord Clamp Cutter



FY2004
Amniotomy Hook



FY2305 Paper/Paper Gusseted Pouch

Size: 90×50×150mm, 90×50×250mm
110×30×190mm, 140×75×250mm etc.



FY2304 Heat-Sealing Gusseted Pouch

Size: (75+25)mm×300mm, (100+50)mm×300mm
(150+50)mm×400mm, (200+55)mm×400mm
(250+60)mm×450mm



FY3903
Surgical Skin Marker

Tip Size: 0.5mm, 1.0mm, 3.5mm
Available with ruler, Sticker
Ink Color: Purple for not easy to wipe out. Blue and Red for easy to wipe out
Gamma Sterile
Can be removed by medical alcohol or low concentration alcohol.



FY2303 Heat-sealing Flat Roll

Size: 55mm×200m, 75mm×200m
100mm×200m, 150mm×200m
400mm×200m



FY2307 Medical Wrapping Paper

Material: Paper or SMS
Color: Blue, Green
Size: 20cm-120cm regular




FY2306 Sterilization Indicator Tape

Size: 12mmx50m, 19mmx50m, 25mmx50m
 Plasma: Turn blue to red under Plasma Sterilization
 Steam: Turn beige to grayish black under Steam Sterilization
 EO : Turn red to green under EO sterilization



FY2309 Indicator Strip

Size: 15mmx100mm, 15mmx200mm
 Plasma: Turn Green to yellow under Plasma Sterilization
 Steam: Class4: Turn blue to Black under Steam Sterilization
 Class5: Turn light blue to dark blue under Steam Sterilization
 Class4: Turn beige to black under Steam Sterilization
 EO: Turn pink to blue under EO sterilization



FY2207 Electro-surgical Pad(with/without Cable)

Size: Adult, Child, Infant
 Type: Mono-polar upright, Mono-polar horizontal
 Bipolar upright, Bipolar horizontal
 Backing: Non-woven or PE Foam



FY2206 Disposable E.C.G Accessories

Material: non-woven, foam, etc
 Size: Adult, Child



FY2210 E.C.G Paper

Classification: single-lead electrocardiograph,
 three-lead electrocardiograph,
 six-lead electrocardiograph,
 twelve-lead electrocardiograph



FY3008 Disposable Oxi-Sensor

9 pin, Neonate/Adult
 Packing: 1pc/blister packing, 40pcs/box. 4boxes/ctn



FY2318 Bowie-Dick Test Pack/Paper

Cat. NO.	Description	Size
FY231801	Bowie-Dick Test Pack	115mmx125mm
FY231802	Bowie-Dick Test Paper	210mmx300mm

Bowie-Dick Pack/paper are pre- assembled, single use test packs designed to detection of defects in the sterilization process. An indicator sheet printed with water based, non-toxic chemical indicator. The color change from yellow to black of indicator after steam sterilization indicates that all air was removed and replace by steam.



FY5001 Emesis Bag

Color: blue/white and ect.
 Material: LDPE
 Tuickness: 30-50micron
 Usage: Used in a variety of situations vomit handing

FY8101 Needle Counter

It can hold needles in place for secure closure and disposal. Foam block needle counters with magnet,
 Size: 10count, 15count, 20count, 40count
 With or without blade remover is available
 Sterile



FY2007 Surgical Light Handle Cover

Disposable foam light handle covers act as a sterile interface between the non-sterile surgical light and the surgeon. The covers provide a simple and economical alternative to re-processing reusable light handles.
 Color: Blue
 EO Sterile



FY5002 Pill Case

Material: High quality PP material
 Features: Waterproof design transparent cover, built-in items at a glance, accurate access, removable grid use and fast.





**FY4001
Ultrasound Gel**

Ultrasound Gel is water soluble, non-greasy and will not irritate tissues or injure the ultrasound applicator or electrodes. It spreads easily and evenly and is an efficient coupling medium for ultrasound transmission. Better sound transmission, distinctive image display, no toxic and irritant effect, no damage and corrosion to the ultrasound probe.
Chemical Composition: Polymer, humectants, water, approved food color, perfume and preservatives.
Size: 60ml, 100ml, 150ml, 250ml



**FY4101
Conductive Gel**

Size: 60ml, 100ml, 150ml, 250ml



FY0619 Tourniquet

Easy to wipe clean
Highly flexible, latex free cotton ribbon which is gentle on the skin.
sturdy snap-on clip
quick and safe releasing



FY0621 Non-latex Tourniquet

Material: TPE



FY0625 Pre-cut Roll Tourniquet

Material: TPE



FY2103 I.D. Bracelet


Material: PVC
Model: write-on, insert-card
Size: infant, adult, adult/infant
Color: blue, pink, white and so on

FY1801 Sharps Container



FY180101

Capacity: 0.8L
Product Size:
Top Dia: 125mm
Bottom Dia: 95mm
Height: 115mm
Carton Size: 510x390x570mm
Quantity in one Carton: 180pcs



FY180102

Capacity: 3.0L
Product Size:
Top Dia: 170mm
Bottom Dia: 150mm
Height: 150mm
Carton Size: 520x520x460mm
Quantity in one Carton: 90pcs



FY180103

Capacity: 3.0L
Product Size:
Top Dia: 200mm
Bottom Dia: 160mm
Height: 140mm
Carton Size: 620x340x450mm
Quantity in one Carton: 70pcs



FY180104

Capacity: 2.0L
Product Size:
Top Dia: 200mm
Bottom Dia: 120mm
Height: 170mm
Carton Size: 405x270x450mm
Quantity in one Carton: 30pcs



FY180105

Capacity: 1.0L
Product Size:
Top Dia: 110mm
Bottom Dia: 90mm
Height: 160mm
Carton Size: 360x320x460mm
Quantity in one Carton: 100pcs



FY180108

Capacity: 3.0L
Product Size:
Top Dia: 150mm
Bottom Dia: 120mm
Height: 260mm
Carton Size: 480x320x480mm
Quantity in one Carton: 30pcs



**FY180112
Needle Cutter**

Item: Needle Cutter
Material: ABS + Blade
Description: Dia 11.33cm x 9.6cm, manual
Packing: 60pcs/ctn, 48x36x49cm



**FY180109
Paper Safety Box**

Size: 2.5L, 5L, etc.



**FY1802
Biohazard Specimen Bag**

Size: 8"x10"
Thickness : 2mil 3ply
Size: 6"x10"
Thickness: 2mil 3ply