

Foyomed[®]

NINGBO FOYOMED MEDICAL INSTRUMENTS CO., LTD.

ADD: 9F, Ningshing Building, No.666, Jiangnan Road,
Hi-Tech Zone, Ningbo, China
Tel: 86-574-81877253
<http://www.foyomed.com> www.foyocare.com
E-mail: products@foyomed.com

MFLAB[®]

NINGBO MFLAB MEDICAL INSTRUMENTS CO., LTD.

ADD: NO. 508, Yindong Road (N) , Yinzhou Economic Development zone
Ningbo, 315145, Zhejiang P.R. China
Tel: 86-574-81877852
<http://www.breathcare.cn> www.mflab.com
E-mail: info@mflab.com



SILICONE MEDICAL PRODUCTS

NINGBO FOYOMED MEDICAL INSTRUMENTS

Establishment

FOYOMED is established in the year 2003 as a professional import and export corporation for medical products. It is the headquarters of MFLAB Medical with production area covering 15,500 m2. Among it, it is including 100000 grade purification workshops and cleaning room, totally 9,200 m2. There are more than 250 workers in the production line.

Team Work

FOYOMED has a high efficient and professional team, full of wisdom and passion, with 15 years of working experience. We take the diligent, professional and responsible spirit, committed to provide customers good products and service.

Export Reputation

Our FOYOMED products have been exported to over 100 countries worldwide through a network of distributors and subsidiaries, covering North American, South American, European, Asian and African countries.

Principle

We FOYOMED adheres to "For Your Smile" principle to strive to achieve excellent quality, responsible price and first-rate service in operation, and cherish every business opportunity.

Leading Products

FOYOMED has a wide range of products in business with leading products respiratory products, disposable medical products, medical equipment and laboratory products.

Brands

There are five brands for our FOYOMED: Foyomed, Foyocare, FoyoLab, Mflab and Smilecare in distribution.

Quality System

FOYOMED is the headauarters of MELAB Medical audited by TUV Rheinland every year, more than 30 products has got CE certificates issued by TUV Rheinland, and also passed by GMP, ISO13485, US FDA on-site Audit, Saudi Arabian FDA on-site Audit, Malaysian MOH on-site Audit, and Korea MOH on-site Audit. We have also got a number of patents.

LB3110 Silicone Manual Resuscitator

- The set with oxygen reservoir bag, Silicone mask and oxygen tubing
- 100% medical level silicone Material, Latex free
- Additional accessories (PEEP Valve, B Valve, etc.) and customized label are available
- All parts(excluding the oxygen tubing & oxygen reservoir) are fully autoclavable up to 134°C
- Reusable for 20 times
- Size: Adult/Pediatric/Infant



LB3110L

Liquid Silicone Manual Resuscitator (Liquid Silicone has better elasticity)

- ◉ The set with oxygen reservoir bag, Silicone mask and oxygen tubing
- ◉ 100% medical level silicone Material, Latex free
- ◉ Additional accessories (PEEP Valve, B Valve, etc.) and customized label are available
- ◉ All parts(excluding the oxygen tubing & oxygen reservoir) are fully autoclavable up to 134°C
- ◉ Reusable for 20 times
- ◉ Size: Adult/Pediatric/Infant



LB3030

Anesthesia Mask (one piece silicone mask)

- ◉ All-in-one piece soft silicone cushion for better seal and comfort
- ◉ Ergonomic shape for lower dead space and better air circulation
- ◉ 100% Medical level Silicone material, Latex free
- ◉ Reusable all parts are fully autoclavable up to 134°C
- ◉ Available in 7 sizes

Cat.NO.	Size	Connector (ID.)	Package
LB303000	0#, Neonate	15mm OD	Ploy bag
LB303001	1#, Infant	15mm OD	Ploy bag
LB303002	2#, Child	22mm ID	Ploy bag
LB303003	3#, Adult S	22mm ID	Ploy bag
LB303004	4#, Adult M	22mm ID	Ploy bag
LB303005	5#, Adult L	22mm ID	Ploy bag
LB303006	6#, Adult XL	22mm ID	Ploy bag



LB3031

Anesthesia Mask (rounded silicone mask)

- ◉ All-in-one piece soft silicone cushion for better seal and comfort
- ◉ Ergonomic shape for lower dead space and better air circulation
- ◉ Made of 100% clear, medical grade silicone, latex free
- ◉ Reusable, all parts are fully autoclavable up to 134°C
- ◉ Packing: Bulk or PE pack available

Cat.NO.	Size	Connector (ID.)	Package
LB303130	000#, Premature	15mm OD	Ploy bag
LB303120	00#, Premature	15mm OD	Ploy bag
LB303110	0#, Neonate	15mm OD	Ploy bag
LB303101	1#, Infant	15mm OD	Ploy bag
LB303102	2#, Child	22mm ID	Ploy bag





**LB3032
Anesthesia Mask (two piece silicone mask)**

- Soft silicone cushion encased in hard plastic shell for ease to handle and better seal
- 100% Medical level Silicone and PC material, Latex free
- Reusable all parts are fully autoclavable up to 134°C
- Available in 3 sizes

Cat.NO.	Size	Connector (ID.)	Package
LB303203	3#,Adult S	22mm	Ploy bag
LB303204	4#,Adult M	22mm	Ploy bag
LB303205	5#,Adult L	22mm	Ploy bag



**LB3033
Manual Sputum Processor**

Cat.NO.	Size	Package
LB303301	0#	Ploy bag
LB303302	1#	Ploy bag
LB303303	2#	Ploy bag

**LB3210
Disposable Silicone Laryngeal Mask**

- Medical grade silicone for super compatibility and patient comfort
- Available with seven sizes for all patient use

Cat.NO.	Type	Applicable People
LB321010	1.0	< 5kg
LB321015	1.5	5-10kg
LB321020	2.0	10-20kg
LB321025	2.5	20-30kg
LB321030	3.0	30-50kg
LB321040	4.0	50-70kg
LB321050	5.0	> 70kg



**LB3220
Reusable Silicone Laryngeal Mask**

- Reusable Silicone Laryngeal Mask are designed to be sterilized and used up to 40 times
- Durable high temperature resistant material suitable for autoclave sterilization
- Additional user card for record of device usage times and patient information is supplied in each pouch


Cat.NO.	Type	Applicable People
LB322010	1.0	< 5kg
LB322015	1.5	5-10kg
LB322020	2.0	10-20kg
LB322025	2.5	20-30kg
LB322030	3.0	30-50kg
LB322040	4.0	50-70kg
LB322050	5.0	>70kg



**LB3212R
Disposable Double Lumen Reinforced Silicone Laryngeal Mask**

- Upgrade Double-Lumen Laryngeal Mask Airway
- Maintaining an effective airway
- Drainage Lumen for suction the Gastric juice suction
- Provides more airway security

Cat.NO.	Type	Applicable People
LB321210R	1.0	< 5kg
LB321215R	1.5	5-10kg
LB321220R	2.0	10-20kg
LB321225R	2.5	20-30kg
LB321230R	3.0	30-50kg
LB321240R	4.0	50-70kg
LB321250R	5.0	> 70kg




LB3213 Disposable Curved Double Lumen Silicone LMA

Product Features

- ◉ Made of imported medical grade silicone, the laryngeal mask is with very good biocompatibility.
- ◉ With high volume cuff design, the cuff of the laryngeal mask fits with patient's oropharyngeal structure properly, which reduces the pressure to patient's pharyngeal mucosa, and improves comfort.
- ◉ The C-shape design based on the physiological bending of human throat makes the insertion easier, faster and more accurate.
- ◉ The laryngeal mask is equipped with gastric drainage channel. Gastric tube can be inserted through it so as to reduce the risk of reflux and aspiration.
- ◉ The cuff is integrated and its sealing pressure is over 35cm H₂O.

Size	Patient Weight	Maximum Size of Gastric Tube	Maximum Volume of Cuff	Maximum Pressure inside Cuff
1.0	0-5KG	6Fr	5ml	≤ 60cmH ₂ O
1.5	5-10KG	8Fr	9ml	≤ 60cmH ₂ O
2.0	10-20KG	10Fr	12ml	≤ 60cmH ₂ O
2.5	20-30KG	12Fr	20ml	≤ 60cmH ₂ O
3.0	30-50KG	14Fr	30ml	≤ 60cmH ₂ O
4.0	50-70KG	14Fr	40ml	≤ 60cmH ₂ O
5.0	70-100KG	14Fr	50ml	≤ 60cmH ₂ O



LB5040S Silicone Reinforced Endotracheal Tube(Cuffed)

Cat. No.	I.D.(mm)	O.D.(mm)
LB504035S	3.5	5.3
LB504040S	4.0	6.0
LB504045S	4.5	6.7
LB504050S	5.0	7.3
LB504055S	5.5	8.0
LB504060S	6.0	8.7
LB504065S	6.5	9.3
LB504070S	7.0	10.0
LB504075S	7.5	10.7
LB504080S	8.0	11.3
LB504085S	8.5	12.0
LB504090S	9.0	12.7
LB504095S	9.5	13.3



Made of medical grade silicone, better biocompatibility
 High pressure silicone cuff maintains good sealing
 Silicone tube with spiral reinforced wire reducing risk of kinking
 100% medical grade silicone and soft tip, reduce the damage to the tracheal wall and mucosa
 Single-use and reusable 2 types available, both can be equipped with a stylet
 Reusable type is suitable for autoclave sterilization only, can be reused up to 10 times


LB3410S Silicone Nasopharyngeal Airway

Features:

- Rounded, beveled tip-promotes safe insertion and improves patient comfort
- Large trumpet provides secure placement
- FR marked on the airway for easy size selection
- Disposable
- Individual packaging
- Latex free
- Sterile

Specifications:

- Trumpet Type, I.D. 3.0-9.0
- Adjustable Type, I.D. 2.5-8.5



LB432S Breathing Circuit for Ventilator


Made of high-quality imported silicone materials and with the method of platinum catalytic technology.

Patent technology of integrated molding, high strength, not prone to coming off and being separated, good flexibility.

Can endure high temperature and high pressure disinfection, reusable.

Joints are made with injection molding with no gas leakage.

Applicable to all kinds of imported and domestic breathing machine, anesthesia machine. Length could be customized and free options among glass, Y type joint, L-shaped connector, masks, breathing air.



**LB431S
Breathing Tube for Anesthesia Apparatus**

Made of high-quality imported silicone materials and with the method of platinum catalytic technology.

Patent technology of integrated molding, high strength, not prone to coming off and being separated, good flexibility.

Can endure high temperature and high pressure disinfection, reusable.

Joints are made with injection molding with no gas leakage.

Applicable to all kinds of imported and domestic breathing machine, anesthesia machine. Length could be customized and free options among glass, Y type joint, L-shaped connector, masks, breathing air.

**LB4930
Anesthesia Reservoir - Silicone**

Cat.NO.	Size
LB493000	0.5L
LB493001	1L
LB493002	2L
LB493003	3L





**LB3041
CPAP Mask**

- ⊙ Size: S, M, L
- ⊙ Made of high-quality liquid medical silicone and suitable for facial features. It will not be aging, hardening but comfortable and durable
- ⊙ The contoured sealing edge head belt eliminate discomfort caused by the oppression to facial nerve
- ⊙ The design of nasal part fits well with the nose and leaves little air blind angle, which limits breath of exhaust and improve comfort
- ⊙ Comfortable, durable silicone skirt design, five points head belt, anti-fog design of standard nose and mouth cover
- ⊙ The sizes which can meet almost all type of faces and do not leak

**LB6220S
Silicone Suction Catheter**



- ⊙ Size: Fr6-Fr20
- ⊙ Made of high-quality imported silicone rubber materials and with the method of platinum catalytic technology.
- ⊙ The tube is soft and without stimulation. There is smooth plug at the top to avoid injury caused respiratory. Damage and bleeding, significantly reducing patients' pain.
- ⊙ Obvious scale marks, easy to know the depth of insertion.
- ⊙ Applicable to sputum suction in various clinical departments.



**LB6410S
Silicone Stomach Tube**

- ⊙ Made of 100% medical grade silicone, the tube is soft and clear.
- ⊙ X-ray opaque line through the whole tube.
- ⊙ Clear depth marks on the tube.
- ⊙ Size: Fr6-Fr18



**LB6420S
Silicone Ryle's Stomach Tube**

- ⊙ Made of 100% medical grade silicone, the tube is soft and clear.
- ⊙ X-ray opaque line through the whole tube.
- ⊙ Clear depth marks on the tube.
- ⊙ Stainless steel balls are sealed in the tip end; it is easier to indwell the tube at a certain position.
- ⊙ Size: Fr6-Fr18



**LB6123
Mucus Specimen Container**

Cat. NO.	Specification	Quantity (pcs)	Measurement(cm)
LB612340	40cc	50	40*27*20
LB612380	80cc	50	40*27*20

For sputum specimen collection, can form closed sputum specimen collection and transportation, clean, convenient, sanitary

All Silicone Foley Catheter

- 100% silicone
- Color-coded for easy size identification
- Alternative types are available
- Sterile by EO
- For single use only
- Plastic valves only

FY0105 2-Way All Silicone Foley Catheter

Cat. NO.	Size (Fr/Ch)	Balloon Capacity (cc)	Color Code
FY010506	6	3	Light Red
FY010508	8	3-5	Black
FY010510	10	3-5	Grey
FY010512	12	5-10/30	White
FY010514	14	5-10/30	Green
FY010516	16	30	Orange
FY010518	18	30	Red
FY010520	20	30	Yellow
FY010522	22	30	Violet
FY010524	24	30	Blue
FY010526	26	30	Pink




All Silicone Foley Catheter

- 100% silicone
- Color-coded for easy size identification
- Alternative types are available
- Sterile by EO
- For single use only
- Plastic valves only

FY0106 3-Way All Silicone Foley Catheter

Cat. NO.	Size (Fr/Ch)	Balloon Capacity (cc)	Color Code
FY010616	16	30	Orange
FY010618	18	30	Red
FY010620	20	30	Yellow
FY010622	22	30	Violet
FY010624	24	30	Blue
FY010626	26	30	Pink





FY0115 Silicone Haematuria Catheter

Features:

- Made of 100% medical grade silicone, safer and less allergy to human body
- Special "anti-blockage" technology enables barrier-free catheterization
- Central positioned high performance balloon
- Providing continuous irrigation, lagger drainage eyes, faster irrigation flow
- Specially designed for the drainage of blood clots & tissues after prostate surgery
- More rigid catheter shaft, better anti-collapsing than latex type
- Equipped with irrigation cap to prevent the urine backflow
- Available with different tips, Straight, Couvelaire, Dufour, etc
- Transparent catheter allow clear observation of irrigation effect changes

Irrigation Time Required For 3000ml Normal Saline




2-Way Standard

Cat.No.				Size (Fr/Ch)	Balloon Volume (ml/cc)	Color Code	Length (mm)
Straight	Couvelaire	Dufour	Mercier				
FY0115216S	FY0115216C	FY0115216D	FY0115216M	16	5-15	Orange	400
FY0115218S	FY0115218C	FY0115218D	FY0115218M	18	50	Red	400
FY0115220S	FY0115220C	FY0115220D	FY0115220M	20	50	Yellow	400
FY0115222S	FY0115222C	FY0115222D	FY0115222M	22	50	Purple	400
FY0115224S	FY0115224C	FY0115224D	FY0115224M	24	50	Blue	400

3-Way Standard

Cat.No.				Size (Fr/Ch)	Balloon Volume (ml/cc)	Color Code	Length (mm)
Straight	Couvelaire	Dufour	Mercier				
FY0115316S	FY0115316C	FY0115316D	FY0115316M	16	5-15	Orange	400
FY0115318S	FY0115318C	FY0115318D	FY0115318M	18	50	Red	400
FY0115320S	FY0115320C	FY0115320D	FY0115320M	20	50	Yellow	400
FY0115322S	FY0115322C	FY0115322D	FY0115322M	22	50	Purple	400
FY0115324S	FY0115324C	FY0115324D	FY0115324M	24	50	Blue	400



FY0122 All Silicone Foley Catheter with Temperature Sensor

- All silicone foley catheter with temperature sensor is designed to be used before, during, and after surgery to monitor core body temperature.
- It is made of 100% medical grade silicone, and is available with long term placement.
- We also offer a various kind of temperature extension cable for different indication.

Cat. NO.	Size (Fr/Ch)	Balloon Capacity (cc)	Color Code
FY012208	8	5	Black
FY012210	10	5	Grey
FY012212	12	10	White
FY012214	14	10	Green
FY012216	16	30	Orange
FY012218	18	30	Red
FY012220	20	30	Yellow
FY012222	22	30	Violet
FY012224	24	30	Blue
FY012226	26	30	Pink

FY0111S
Non-Adhesive Silicone Male External Catheter

- ⊙ Soft, breathable silicone is designed to be non-irritating and can be crimped for a custom fit
- ⊙ Transparent material allows for skin health monitoring while the catheter is in place
- ⊙ Skin-friendly adhesive tape helps prevent accidental detachment
- ⊙ Not made with natural rubber latex


Cat.No.	Type	Diameter				
		S	M	L	XL	XXL
FY0111S	Standard	24mm	28mm	31mm	35mm	40mm



FY0121
Self Adhesive Silicone Male External Catheter

- ⊙ Made of 100% imported medical grade silicone
- ⊙ Pop On: Suitable for old person or sports person
- ⊙ Standard: Suitable for regular length and adhesive for daily wear
- ⊙ Wide Band: More adhesive and viscous force than standard for added security

Type	Diameter				
	S	M	L	XL	XXL
Pop On	24mm	28mm	31mm	35mm	40mm
Standard	24mm	28mm	31mm	35mm	40mm
Wide Band	24mm	28mm	31mm	35mm	40mm




**LB6125
Silicone Gastrostomy Tube**

Gastrostomy tube is usually used for patients with advanced esophageal cancer. Doctors do not remove the cancer, but only leave a tube in the abdomen. The purpose of the gastrostomy tube is to provide a long-term enteral nutrition for patients who have normal stomach function but cannot eat food. This treatment is relatively simple and safe and economical.

- ⊙ Made of 100% medical grade silicone, the tube is soft and clear
- ⊙ X-Ray opaque line through the whole tube
- ⊙ The balloon is glued to the main tube both inside and outside, it is elastic and flexible.


Cat. No.	Size(Fr)	Color Code	Length (mm)	Marked Position(cm)
LB612512	12Fr/5ml	White	240	2/4/6
LB612514	14Fr/5ml	Green	240	2/4/6
LB612516	16Fr/5ml	Orange	240	2/4/6
LB612518	18Fr/20ml	Red	240	2/4/6
LB612520	20Fr/20ml	Yellow	240	2/4/6
LB612522	22Fr/20ml	Purple	240	2/4/6
LB612524	24Fr/20ml	Dark Blue	240	2/4/6




**LB7010S
Silicone Chest Drainage Tube**

- ⊙ Made of 100% medical grade silicone
- ⊙ X-ray opaque line through the whole tube
- ⊙ Size: Fr12-Fr36

FY0109S
Silicone Penrose Drainage Tube




- ⊙ 100% medical grade silicone for superior biocompatibility
- ⊙ X-Ray opaque line throughout the length
- ⊙ Soft and smooth tube minimizes trauma, tooth-like structure to minimize clog or obstruct
- ⊙ Clear connecting tube provides insoection and drainage visibility
- ⊙ Sterilized double package
- ⊙ Size: 6mm x 30cm, 6mm x 50cm, 8mm x 30cm, 8mm x 50cm, 10mm x 30cm, 10mm x 50cm, 12mm x 30cm, 12mm x 50cm



FY0120
T-shaped Silicone Perforated Drain

- ⊙ Made of 100% imported medical grade silicone.
- ⊙ With X-ray Opaque and super biocompatibility, it can be used for wound draining.

Cat. NO.	Size(Fr)	O.D.(mm)
SH025612	12	4.0
SH025614	14	4.7
SH025616	16	5.3



FY0235
Closed Wound Drainage System Kit

- Made of 100% imported medical grade silicone
- This product includes 100cc/200cc silicone reservoir, 7mm/10mm Flat drains or other round drains and stainless steel trocar
- Can use separately or using trocar for insertion and combining with silicone Reservoir for fluids syction and collection
- Used for various surgeries
- Supplied sterile, EO sterilization



FY0236
Silicone Reservoir/Bulb

- Made of 100% imported medical-grade silicone, with anti-reflyx valve.
- Can compatible with any of following drains: Flat Drain, Fluted Drain and Round Drain.
- Squeeze the reservoir and connect with drains for suction, the low-vacuum reservoir can make stable and effective suction for exudates.

Cat. NO.	Size(cc)	Marked Position(ml)
FY023601	100	25
FY023602	200	30
FY023603	400	40

*Clamp available.


FY0243
Silicone Round Channel/Fluted Drain

- Made of 100% imported medical grade silicone.
- With X-ray opaque, used for closed wound suction draining.
- These intrnal four channels, or fluted, are utilized to quickly remove exudates.
- Using one-time extrusion, and hubless of the tube minimizes pain associated with device removal.

Cat. NO.	Size(Fr)	O.D.(mm)	Length(cm)
FY024310	10	3.3	120
FY024315	15	5.0	120
FY024319	19	6.3	120
FY024324	24	8.0	120

*Length could be customized.





**FY0244
Silicone Round Perforated Drain**

- Made of 100% imported medical grade silicone.
- X-ray opaque line through the whole tube.
- It can be equipped with reservoir and trocar.
- Clear tubing provides, inspection and drainage visibility.

Cat. NO.	Size(Fr)	O.D.(mm)
FY024410	10	3.0
FY024412	12	4.0
FY024414	14	4.7
FY024416	16	5.3
FY024418	18	6.0
FY024420	20	6.7

*Length could be customized.


**FY0245
Silicone Flat Perforated Drain**

- Made of 100% imported medical grade silicone.
- This Drian combines a transparent silicone round tube and a white flat drain.
- The inner wall of the implanted flat segment features tooth-like to prevent possible occlusion when tube is pressed or kinked.
- 3/4 or full perforation are available for optimal drainage with X-ray opaque effectiveness throughout its entire implanted segment.
- Using Trocar for insertion and combining with silicone. Reservoir for fluids suction and collection.

Cat. NO.	Width(mm)	Length(cm)
FY024507	7	25
FY024510	10	30
FY024513	13	34

*Length could be customized.





**FY0242
Silicone Flat Fluted Drain**

- Made of 100% imported medical grade silicone.
- This Drain combines a transparent silicone round tube and a white flat drain. The inner wall of the implanted flat segment features multiple narrow ducts to prevent possible occlusion when the tube is pressed or kinked.
- Full fluted design is available for optimal drainage with X-ray opaque effectiveness throughout its entire implanted segment.
- Using Trocar for insertion and combining with silicone reservoir for fluids suction and collection.

Cat. NO.	Width(mm)	Length(cm)
FY024208	8	30
FY024211	11	30

*Length could be customized.



Catheter (Fr)	Length (cm)	Balloon Volume (ml)
18	40	80

**FY0323
Disposable Cervical Ripening Balloon**

Introduce:

The Cervical Ripening Balloon is used for physical dilation of cervical canal prior while the cervix is unfavorable for induction. Balloons are inflated with sterilized saline.

By using the pressure of both balloon, the cervix will synthesize and release prostaglandin, and promote the cervix ripening.

Features:

- Shorten the process of birth, reduce the pain of parturition
- Medical grade material ensures the highest level of safety
- Obvious effect on promoting the cervix ripening
- Ripens and dilates the cervix without pharmaceuticals
- Silicone body material eliminates the issue of patient sensitively to latex



Catheter (Fr)	Length (cm)	Balloon Volume (ml)
24	55	500

FY032401 Postpartum Balloon

Introduce:

Used for Temporary control or reduction of postpartum uterine bleeding when conservative management is warranted.

For while the balloon was inflated with sterilized saline, the postpartum balloon applies pressure to the uterine walls that will help to hemostasis.

Features:

- ⊙ Easy to operate and position
- ⊙ Safe and efficient way to promote hemostasis
- ⊙ Can be used for both vaginal and cesarean delivery
- ⊙ Silicone body material eliminates the issue of patient sensitively to latex

FY032403 Balloon Uterine Stent

Introduce:

- ⊙ Used for placement in the uterus following an intrauterine surgical procedure in order to prevent intrauterine adhesious and reduce uterine bleeding
- ⊙ Obvious effect on intrauterine adhesious prevention
- ⊙ Excellent way to reduce uterine bleeding
- ⊙ No contraindication and after effect


Cat.No.	Catheter(Fr)	Balloon Length(cm)	Balloon Width(cm)	Max.Volume(ml)
FY03240301	9	3	2.8	5
FY03240302	9	4	4	8



FY0322 HSG Catheter

Introduce:

- Medical silicone material, good biocompatibility.
- The balloon wall is thin yet strong and smooth, while the connection between the balloon and the tube body is also smooth. These features effectively reduce the risk of trauma to the lumen during insertion and removal, thereby enhancing patient comfort.
- The catheter body is designed with "X-rays" to allow for easy observation of the catheter under imaging equipment within the body.
- Additionally, it comes with a guide wire that reinforces the catheter's hardness, making it easier to introduce into the fallopian tube lumen.



Catheter(Fr)	Length(mm)	Balloon Volume(ml)
5	360	1.5
7	360	1.5



FY0628 Silicone Loop-Vessel Loops

Features:

- Identification: red for artery, blue for vein, yellow for nerve, and white for tendon
- Made from durable medical grade silicone
- Fluid-repellent and neatly wrapped to prevent tangling
- Radiopaque, sterilized, latex free

Micro

Cat.No.	Width(mm)	Thickness(mm)	Length(mm)	Colors
FY062801B	0.7	0.5	460	Blue
FY062801R	0.7	0.5	460	Red
FY062801W	0.7	0.5	460	White
FY062801Y	0.7	0.5	460	Yellow

Mini

Cat.No.	Width(mm)	Thickness(mm)	Length(mm)	Colors
FY062802B	1.33	1.05	460	Blue
FY062802R	1.33	1.05	460	Red
FY062802W	1.33	1.05	460	White
FY062802Y	1.33	1.05	460	Yellow

Maxi

Cat.No.	Width(mm)	Thickness(mm)	Length(mm)	Colors
FY062803B	2.4	1.22	460	Blue
FY062803R	2.4	1.22	460	Red
FY062803W	2.4	1.22	460	White
FY062803Y	2.4	1.22	460	Yellow

Super Maxi

Cat.No.	Width(mm)	Thickness(mm)	Length(mm)	Colors
FY062804B	5.3	1.5	460	Blue
FY062804R	5.3	1.5	460	Red
FY062804W	5.3	1.5	460	White
FY062804Y	5.3	1.5	460	Yellow


FY13253 Silicone Foam Dressing with Border

Product Features:

- ⊙ Gentle on the skin and painless to remove, with no residue, waterproof, soft and comfortable.
- ⊙ "Vertically" absorbs wound mucosal exudate, reducing the risk of maceration around the wound, absorption capacity is great.
- ⊙ Can be reused after inspection of the wound and will not adhere to the wound bed.
- ⊙ Helps to reduce scarring.

Size: 3"x3", 4"x4", 5"x5", 6"x6", 6"x8"

Other size could be customized



Water-proof PU Film

Silicone layer with holes

Hyper absorbent foam

Isolation film

FY1222 Silicone Tape

Product Features:

- ⊙ Gentle on the skin and painless to remove.
- ⊙ Reusable, does not stick to hair and leaves no residue, no remains.
- ⊙ Easy to use and hypoallergenic, with good adhesion.
- ⊙ Waterproof, but permeable.

Intended Use:

- ⊙ Terrific for use on patients where skin condition is fragile, poor or damaged, such as patients after radiation therapy, chemotherapy, in dialysis or infants.
- ⊙ Applicable to a variety of in-dwelling catheters, electrodes, or other external fixation that requires gentle, but robust skin protection and adhesion.



FY0310 Menstrual Cup

This menstrual cup is reusable, made of medical grade silicone, which is widely applied in medical care products such as infant nipples, medical catheters etc, because it is toxic free and allergy free, safe to use. It is designed to be folded inserted into your vagina to collect your menstrual blood. Once inside, it opens up and sits lower than tampon. It is held in place by your vagina muscles.

Advantages:

Keep cool and safe

Comfortable, clean and easy to use

100% medical grade silicone, no BPA or latex

Reusable, eco-friendly and economical

Leak-free protection for up to 10 hours at a time

Long-term use can reduce the risk of gynecological inflammation

Worry-free when traveling, swimming or exercising during menstrual period

S: Recommended for women who have never delivered vaginally or have delivered by caesarean section

L: Recommended for women who have delivered vaginally

